The Wound Healing Company

SUPRATHEL® 250

INSTRUCTIONS FOR USE

- EN Alloplastic skin substitute
- DE Alloplastischer Hautersatz
- BG Алопластичен кожен заместител
- CS Alopplastická kožní náhrada
- DA Alloplastisk huderstatning
- EL Αλλοπλαστικό υποκατάστατο δέρματος
- ES Reemplazo aloplástico de la piel
- ET Alloplastikas kasutatav naha aseaine
- FI Alloplastinen ihon korvike
- FR Substitut alloplastique de la peau
- GA Leaschraiceann allaplaisteach
- HR Aloplastični nadomjestak za kožu
- HU Aloplastikus bőrhelyettesítő
- IT Sostituto alloplastico della cute
- LT Aloplastinis odos pakaitalas
- LV Aloplastisks ādas aizvietotājs
- MT Sostitut tal-ġilda alloplastiku
- NL Alloplastische huidvervanger
- NO Alloplastisk huderstatning
- PL Alloplastycznym substytut skóry
- PT Substituto aloplástico da pele
- RO Inlocuitor aloplastic al pielii
- SK Alopplastická náhrada kože
- SL Aloplastični nadomestek kože
- SV Alloplastiskt hud substitut
- TR Alloplastik deri ikamesi

MANUFACTURED BY
PolyMedics Innovations GmbH

EN Instructions for Use SUPRATHEL®250

Product description

Suprathel®250 is an absorbable, microporous membrane and approved as alloplastic skin substitute for the treatment of epidermal and dermal wounds.

Composition

Copolymer composed out of lactide (73–80 %), trimethylene carbonate (12–16 %) and caprolactone (8–11 %).

Indications

Suprathel®250 is used for patients with epidermal and dermal wounds, such as abrasions, split skin graft donor sites, 2nd degree burns as well as 2nd degree burns mixed with 3rd degree burned areas. Suprathel®250 is used for patients with chronic wounds, such as venous and arterial ulcers, as well as diabetic wounds.

Contraindications

Suprathel®250 should not be used on infected wound sites or on severe bleeding wounds without additional hemostatic treatment. Suprathel®250 should not be applied on chronic dry wounds.

Properties and Mechanism

Suprathel®250 is an alloplastic, absorbable skin substitute that is highly permeable to oxygen and water vapour. Suprathel®250 can be adapted to any body part due to its elastic/plastic properties.

Application Procedure

Clean and debride the wound prior to applying Suprathel®250. Remove Suprathel®250 from the package in an aseptic manner. Be sure to separate and discard the green colored protective paper that surrounds the white Suprathel®250 membrane. Apply Suprathel®250 on the affected area in either direction or side. If multiple dressings are required to cover larger areas, Suprathel®250 may be applied overlapping to assure complete coverage of the wound area. Slight stretching of Suprathel®250 may facilitate contouring of the dressing around curved anatomies such as elbows, fingers, etc. Assure complete contact of Suprathel®250 with the affected area by removing all air bubbles and wrinkles. Suprathel®250 may be trimmed so that only injured tissues are covered by the dressing. For protection against mechanical

disruptions, apply a single sheet of fatty gauze or similar dressing material over Suprathel®250.

Aftercare Procedure

A gauze dressing applied over the fatty gauze can act as a protective barrier against mechanical alteration or contamination. This dressing can be changed as needed. Removal of Suprathel®250 during such dressing changes should be avoided to prevent disruption of the underlying healing tissues below the dressing. In case of chronic wounds, depending on the condition of the wound, it may be necessary to change or reapply Suprathel®250 at certain intervals. The transparent nature of applied Suprathel®250 facilitates visualisation of healing tissue.

Adverse reactions

There have been no serious adverse side effects attributed to the product.

Storage and Handling

Suprathel®250 may only be stored in original packaging dry and cool (-20 up to +22 °C). Keep the product away from heat sources. Suprathel®250 can be disposed of with household waste.

Warnings and Precautions

Do not apply a product, where the sterility may not be ensured as this may lead to severe infections. The content is sterile unless sterile packaging is damaged. In case of packaging damages, the sterility of the product is not ensured. The unused contents of opened or damaged sterile packages are to be discarded. Do not reuse and do not resterilise. If the product is nevertheless reused, this may lead to impairment of product performance characteristics (reduced permeability, elasticity, adherence capability as well as sterility). Such changes of material properties may in turn lead to treatment impairments, such as inadequate wound healing as well as infections. In the case of known allergies against components of Suprathel®250, the membrane should not be applied. Suprathel®250 should be removed immediately if there are any signs of allergic reactions to the material. Suprathel®250 should be removed in cases of severe pain or accumulations of wound secretions. Coverage of intact skin may lead to skin macerations and should be avoided.


Additional Information

Suprathel®250 is gamma sterilised. Suprathel®250 is non-toxic. To be applied by healthcare professionals only. Suprathel®250 is a white membrane. In the case of a transparent membrane appearance please do not apply. As Suprathel®250 is a biosynthetic product, slight variations in shape and size may occur. These variations can be compensated by slight manual stretching of the membrane, if required. The summary of safety and clinical performance for Suprathel®250 is publicly available in the European Medical Device Database (Eudamed: <https://ec.europa.eu/tools/eudamed>) via the Basic UDI-DI 426018402AAA0000001PQ or via our homepage (<https://polymedics.com>).

Notice: any serious incident that has occurred in relation to Suprathel®250 should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. Users should inform the patients of residual risks associated with the medical device, particularly with respect to potentially serious infections, as indicated in these instructions for use at section contraindications and warnings and precautions.

This IFU applies to the REF numbers of Suprathel®250, including 110505-250-MDR, 150505-250-MDR, 110910-250-MDR, 150910-250-MDR, 111810-250-MDR, 151810-250-MDR, 111823-250-MDR, 151823-250-MDR.

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PolyMedics Innovations GmbH (PMI)