INSTRUCTIONS FOR USE

DE **Gebrauchsinformation**

EN **Instructions for use**

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TR **Kullanma talimatı**

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SE **Bruksanvisning**

NO **Bruerveiledning**

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DA **Brugerinformation**

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
EN Instructions for use

Product description

Suprathel® is an absorbable, microporous membrane and approved as alloplastic skin substitute for the treatment of epidermal and dermal wounds.

Composition

Copolymer composed out of lactide (73 – 80 %), trimethylene carbonate (12 – 16 %) and caprolactone (8 – 11 %).

Indications

Suprathel® is used for patients with epidermal and dermal wounds, such as abrasions, split skin graft donor sites, 2nd degree burns as well as 2nd degree burns mixed with 3rd degree burned areas.

Contraindications

Suprathel® should not be used on infected wound sites or on severe bleeding wounds without additional hemostatic treatment.

Properties and Mechanism

Suprathel® is an alloplastic, absorbable skin substitute that is highly permeable to oxygen and water vapour. Suprathel® can be adapted to any body part due to its elastic/plastic properties.

Application Procedure

Clean and debride the wound prior to applying Suprathel®. Remove Suprathel® from the package in an aseptic manner. Be sure to separate and discard the green colored protective paper that surrounds the white Suprathel® membrane. Apply Suprathel® on the affected area in either direction or side. If multiple dressings are required to cover larger areas, Suprathel® may be applied overlapping to assure complete coverage of the wound area. Slight stretching of Suprathel® may facilitate contouring of the dressing around curved anatomies such as elbows, fingers, etc. Assure complete contact of Suprathel®

with the affected area by removing all air bubbles and wrinkles. Suprathel® may be trimmed so that only injured tissues are covered by the dressing. For protection against mechanical disruptions, apply a single sheet of fatty gauze or similar dressing material over Suprathel®.

Aftercare Procedure

A gauze dressing applied over the fatty gauze can act as a protective barrier against mechanical alteration or contamination. This dressing can be changed as needed. Removal of Suprathel® during such dressing changes should be avoided to prevent disruption of the underlying healing tissues below the dressing. The transparent nature of applied Suprathel® facilitates visualization of healing tissue.

Adverse reactions

There have been no serious adverse side effects attributed to the product. In rare cases very mild inflammatory response typical for foreign body reaction may be observed. This usually declines during degradation of the membrane.

Storage and Handling

Suprathel® may only be stored in original packaging dry and cool (-20 to +22 °C). Keep the product away from heat sources.

Warnings and Precautions

The content is sterile unless sterile packaging is damaged. In case of packaging damages, the sterility of the product is not ensured. The unused contents of opened or damaged sterile packages are to be discarded. Do not reuse and do not resterilise. If the product is nevertheless reused, this may lead to impairment of product performance characteristics (reduced permeability, elasticity, adherence capability as well as sterility). Such changes of material properties may in turn lead to treatment impair-

ments, such as inadequate wound healing as well as infections. In the case of known allergies against components of Suprathel®, the membrane should not be applied. Suprathel® should be removed immediately if there are any signs of allergic reactions to the material. Suprathel® should be removed in cases of severe pain or accumulations of wound secretions.


Additional Information

Suprathel® is gamma sterilised. Suprathel® is non-toxic. To be applied by healthcare professionals only. Suprathel® is a white membrane. In the case of a transparent membrane appearance please do not apply. As Suprathel® is a biosynthetic product, slight variations in shape and size may occur. These variations can be compensated by slight manual stretching of the membrane, if required.

Notice: any serious incident that has occurred in relation to Suprathel® should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

This IFU applies to the REF numbers of Suprathel®, including 150505, 110910, 150910, 111810, 151810, 111823, 151823 and 12H001.

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PolyMedics Innovations GmbH (PMI)