# SUPRA SDRM<sup>®</sup> 1100

## INSTRUCTIONS FOR USE

EN	Alloplastic skin substitute
DE	Alloplastischer Hautersatz
BG	Алопластичен кожен заместител
CS	Aloplastická kožní náhrada
DA	Alloplastisk huderstatning
EL	Αλλοπλαστικό υποκατάστατο δέρματος
ES	Reemplazo aloplástico de la piel
ET	Alloplastikas kasutatav naha aseaine
FI	Alloplastinen ihon korvike
FR	Substitut alloplastique de la peau
GA	Leaschraiceann allaplaisteach
HR	Aloplastični nadomjestak za kožu
HU	Alloplasztikus bőrhelyettesítő
IT	Sostituto alloplastico della cute
LT	Aloplastinis odos pakaitalas
LV	Aloplastisks ādas aizvietotājs
MT	Sostitut tal-ġilda alloplastiku
NL	Alloplastische huidvervanger
NO	Alloplastisk huderstatning
PL	Alloplastyczny substytut skóry
PT	Substituto aloplástico da pele
RO	Inlocuitor aloplastic al pielii
SK	Aloplastická náhrada kože
SL	Aloplastični nadomestek kože
SV	Alloplastiskt hudsubstitut
TR	Alloplastik deri ikamesi

# EN Instructions for Use SUPRA SDRM®1100

## Product description

SupraSDRM®1100 is an absorbable foam membrane and an alloplastic skin substitute for the treatment of epidermal and dermal wounds.

## Composition

Copolymer composed out of lactide (73 – 80 %), trimethylene carbonate (12 – 16 %) and caprolactone (8 – 11 %).

## Indications

SupraSDRM®1100 is used for patients with epidermal and dermal wounds, such as abrasions, split skin graft donor sites, 2<sup>nd</sup> degree burns as well as 2<sup>nd</sup> degree burns mixed with 3<sup>rd</sup> degree burned areas. SupraSDRM®1100 is used for patients with chronic wounds (excluding dry chronic wounds), such as venous and arterial ulcers, as well as diabetic wounds.

## Contraindications

SupraSDRM®1100 should not be used on infected wound sites or on severe bleeding wounds without additional hemostatic treatment. SupraSDRM®1100 should not be applied on dry chronic wounds.

## Properties and Mechanism

SupraSDRM®1100 is an alloplastic, absorbable skin substitute that is highly permeable to oxygen and water vapour. SupraSDRM®1100 can be adapted to any body part due to its elastic/plastic properties.

## Application Procedure

Clean and debride the wound prior to applying SupraSDRM®1100. Remove SupraSDRM®1100 from the package in an aseptic manner. Be sure to separate and discard the green colored protective paper that surrounds the white SupraSDRM®1100 membrane. Apply SupraSDRM®1100 on the affected area in either direction or side. If multiple dressings are required to cover larger areas, SupraSDRM®1100 may be applied overlapping to assure complete coverage of the wound area. Slight stretching of SupraSDRM®1100 may facilitate contouring of the dressing around curved anatomies such as elbows, fingers, etc. Assure complete contact of SupraSDRM®1100 with the affected area by removing all air bubbles and wrinkles. SupraSDRM®1100 may be trimmed so that only injured tissues are

covered by the dressing. For protection against mechanical disruptions, apply a single sheet of fatty gauze or similar dressing material over SupraSDRM®1100.

## Aftercare Procedure

A gauze dressing applied over the fatty gauze can act as a protective barrier against mechanical alteration or contamination. This dressing can be changed as needed. Removal of SupraSDRM®1100 during such dressing changes should be avoided to prevent disruption of the underlying healing tissues below the dressing. In case of chronic wounds, depending on the condition of the wound, it may be necessary to change or reapply SupraSDRM®1100 in appropriate time intervals depending on individual wound requirements. The transparent nature of applied SupraSDRM®1100 facilitates visualisation of healing tissue.

## Adverse reactions

There have been no serious adverse side effects attributed to the product.

## Storage and Handling

SupraSDRM®1100 may only be stored in original packaging dry and cool (-20 to +25 °C). Keep the product away from heat sources. SupraSDRM®1100 can be disposed of with household waste.

## Warnings and Precautions

Do not apply a product, where the sterility may not be ensured as this may lead to severe infections. The content is sterile unless sterile packaging is damaged. In case of packaging damages, the sterility of the product is not ensured. The unused contents of opened or damaged sterile packages are to be discarded. Do not reuse and do not resterilise. If the product is nevertheless reused, this may lead to impairment of product performance characteristics (reduced permeability, elasticity, adherence capability as well as sterility). Such changes of material properties may in turn lead to treatment impairments, such as inadequate wound healing as well as infections. In the case of known allergies against components of SupraSDRM®1100, the membrane should not be applied. SupraSDRM®1100 should be removed immediately if there are any signs of allergic reactions to the material. SupraSDRM®1100

should be removed in cases of severe pain or accumulations of wound secretions. Coverage of intact skin may lead to skin macerations and should be avoided.


### Additional Information

SupraSDRM®1100 is gamma sterilised. SupraSDRM®1100 is non-toxic. To be applied by healthcare professionals only. SupraSDRM®1100 is a white membrane. In the case of a transparent membrane appearance please do not apply. As SupraSDRM®1100 is a biosynthetic product, slight variations in shape and size may occur. These variations can be compensated by slight manual stretching of the membrane, if required. The summary of safety and clinical performance for SupraSDRM®1100 is publicly available in the European Medical Device Database (Eudamed: <https://ec.europa.eu/tools/eudamed>) via the Basic UDI-DI 426018402AAA0000001PQ or via our homepage (<https://polymedics.com>).

Notice: any serious incident that has occurred in relation to SupraSDRM®1100 should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. Users should inform the patients of residual risks associated with the medical device, particularly with respect to potentially serious infections, as indicated in these instructions for use at section contraindications and warnings and precautions.

This IFU applies to the REF numbers of SupraSDRM®1100, including 31D012-1100-MDR, 31D018-1100-MDR, 31D024-1100-MDR, 310101-1100-MDR, 310202-1100-MDR, 310404-1100-MDR, 310505-1100-MDR, 310909-1100-MDR, 310912-1100-MDR, 311809-1100-MDR, 311818-1100-MDR.

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PolyMedics Innovations GmbH (PMI)