# SUPRA SDRM<sup>®</sup> CW

## INSTRUCTIONS FOR USE

DE	Gebrauchsinformation
EN	Instructions for use
FR	Notice d'utilisation
IT	Informazioni per l'uso
ES	Información sobre el uso
PT	Instruções de Utilização
TR	Kullanma talimatı
PL	Instrukcja Obsługi
HU	Használati utasítás
FI	Käyttöohje
SE	Bruksanvisning
NO	Bruerveiledning
NL	Gebruiksaanwijzing
DA	Brugerinformation

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Tel.: +49 (0)711 719 500-0

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**Product description**

SupraSDRM®CW is an absorbable foam membrane and an alloplastic skin substitute for the treatment of epidermal and dermal wounds.

**Composition**

Copolymer composed out of lactide (73 – 80 %), trimethylene carbonate (12 – 16 %) and caprolactone (8 – 11 %).

**Indications**

SupraSDRM®CW is used for chronic wounds, such as venous and arterial ulcers, as well as for diabetic wounds.

**Contraindications**

SupraSDRM®CW should not be applied to infected areas. SupraSDRM®CW should not be applied to chronic dry wounds. SupraSDRM®CW should not be applied on severe bleeding wounds without additional hemostatic treatment.

**Properties and Mechanism**

SupraSDRM®CW is an alloplastic, absorbable skin substitute that is highly permeable to oxygen and water vapour. SupraSDRM®CW can be adapted to any body part due to its elastic/plastic properties.

**Application Procedure**

Clean and debride the wound prior to applying SupraSDRM®CW. Remove SupraSDRM®CW from the package in an aseptic manner. Be sure to separate and discard the green colored protective paper that surrounds the white SupraSDRM®CW membrane. Apply SupraSDRM®CW on the affected area in either direction or side. If multiple dressings are required to cover larger areas, SupraSDRM®CW may be applied overlapping to assure complete coverage of the wound area. Slight stretching of SupraSDRM®CW may facilitate contouring of the dressing around curved anatomies such as elbows, fingers, etc.

Assure complete contact of SupraSDRM®CW with the affected area by removing all air bubbles and wrinkles. SupraSDRM®CW may be trimmed so that only injured tissues are covered by the dressing. For protection against mechanical disruptions, apply a single sheet of fatty gauze or similar dressing material over SupraSDRM®CW.

**Aftercare Procedure**

A gauze dressing applied over the fatty gauze can act as a protective barrier against mechanical alteration or contamination. This dressing can be changed as needed. The frequency of SupraSDRM®CW dressing changes should be determined based on the extent of wound secretion. The transparent nature of applied SupraSDRM®CW facilitates visualization of healing tissue.

**Adverse reactions**

There have been no serious adverse side effects attributed to the product. In rare cases very mild inflammatory response typical for foreign body reaction may be observed. This usually declines during degradation of the membrane.

**Storage and Handling**

SupraSDRM®CW may only be stored in original packaging dry and cool (-20 to +25 °C). Keep the product away from heat sources.

**Warnings and Precautions**

The content is sterile unless sterile packaging is damaged. In case of packaging damages, the sterility of the product is not ensured. The unused contents of opened or damaged sterile packages are to be discarded. Do not reuse and do not resterilise. If the product is nevertheless reused, this may lead to impairment of product performance characteristics (reduced permeability, elasticity, adherence capability as well as sterility). Such changes of material properties may in turn lead to treatment impairments, such as inadequate wound

healing as well as infections. In the case of known allergies against components of SupraSDRM®CW, the membrane should not be applied. SupraSDRM®CW should be removed immediately if there are any signs of allergic reactions to the material. SupraSDRM®CW should be removed in cases of severe pain or accumulations of wound secretions.


### Additional Information

SupraSDRM®CW is gamma sterilised. SupraSDRM®CW is non-toxic. To be applied by healthcare professionals only. SupraSDRM®CW is a white membrane. In the case of a transparent membrane appearance please do not apply. As SupraSDRM®CW is a biosynthetic product, slight variations in shape and size may occur. These variations can be compensated by slight manual stretching of the membrane, if required.

Notice: any serious incident that has occurred in relation to SupraSDRM®CW should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

This IFU applies to the REF numbers of SupraSDRM®CW, including 41D012, 41D018, 41D024, 410101, 410202, 410404, 410505, 410909, 410912, 411809 and 411818.

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PolyMedics Innovations GmbH (PMI)