
SUPRA SDRM®

INSTRUCTIONS FOR USE

- EN Alloplastic skin substitute
- DE Alloplastischer Hautersatz
- BG Алопластичен кожен заместител
- CS Alopplastická kožní náhrada
- DA Alloplastisk huderstatning
- EL Αλλοπλαστικό υποκατάστατο δέρματος
- ES Reemplazo aloplástico de la piel
- ET Alloplastikas kasutatav naha aseaine
- FI Alloplastinen ihon korvike
- FR Substitut alloplastique de la peau
- GA Leaschraiceann allaplaisteach
- HR Alopplastični nadomjestak za kožu
- HU Alopplastikus bőrhelyettesítő
- IT Sostituto alloplastico della cute
- LT Alopplastinis odos pakaitalas
- LV Alopplastisks ādas aizvietotājs
- MT Sostitut tal-ġilda alloplastiku
- NL Alloplastische huidvervanger
- NO Alloplastisk huderstatning
- PL Alopplastycznym substytut skóry
- PT Substituto aloplástico da pele
- RO Inlocuitor aloplastic al pielii
- SK Alopplastická náhrada kože
- SL Alopplastični nadomestek kože
- SV Alopplastiskt hud substitut
- TR Alopplastik deri ikamesi

Product description

SupraSDRM® is an absorbable foam membrane and an alloplastic skin substitute for the treatment of epidermal and dermal wounds.

Composition

Copolymer composed out of lactide (73–80 %), trimethylene carbonate (12–16 %) and caprolactone (8–11 %).

Indications

SupraSDRM® is used for patients with epidermal and dermal wounds, such as abrasions, split skin graft donor sites, 2nd degree burns as well as 2nd degree burns mixed with 3rd degree burned areas. SupraSDRM® is used for patients with chronic wounds, such as venous and arterial ulcers, as well as diabetic wounds.

Contraindications

SupraSDRM® should not be used on infected wound sites or on severe bleeding wounds without additional hemostatic treatment. SupraSDRM® should not be applied on chronic dry wounds.

Properties and Mechanism

SupraSDRM® is an alloplastic, absorbable skin substitute that is highly permeable to oxygen and water vapour. SupraSDRM® can be adapted to any body part due to its elastic/plastic properties.

Application Procedure

Clean and debride the wound prior to applying SupraSDRM®. Remove SupraSDRM® from the package in an aseptic manner. Be sure to separate and discard the green colored protective paper that surrounds the white SupraSDRM® membrane. Apply SupraSDRM® on the affected area in either direction or side. If multiple dressings are required to cover larger areas, SupraSDRM® may be applied overlapping to assure complete coverage of the wound area. Slight stretching of SupraSDRM® may facilitate contouring of the dressing around curved anatomies such as elbows, fingers, etc. Assure complete contact of SupraSDRM® with the affected area by removing all air bubbles and wrinkles. SupraSDRM® may be trimmed so that only injured tissues are covered by the dressing. For protection against mechanical disruptions, apply a single sheet of fatty gauze or similar dressing material

over SupraSDRM®.

Aftercare Procedure

A gauze dressing applied over the fatty gauze can act as a protective barrier against mechanical alteration or contamination. This dressing can be changed as needed. Removal of SupraSDRM® during such dressing changes should be avoided to prevent disruption of the underlying healing tissues below the dressing. In case of chronic wounds, depending on the condition of the wound, it may be necessary to change or reapply SupraSDRM® at certain intervals. The transparent nature of applied SupraSDRM® facilitates visualisation of healing tissue.

Adverse reactions

There have been no serious adverse side effects attributed to the product.

Storage and Handling

SupraSDRM® may only be stored in original packaging dry and cool (-20 to +25 °C). Keep the product away from heat sources. SupraSDRM® can be disposed of with household waste.

Warnings and Precautions

Do not apply a product, where the sterility may not be ensured as this may lead to severe infections. The content is sterile unless sterile packaging is damaged. In case of packaging damages, the sterility of the product is not ensured. The unused contents of opened or damaged sterile packages are to be discarded. Do not reuse and do not resterilise. If the product is nevertheless reused, this may lead to impairment of product performance characteristics (reduced permeability, elasticity, adherence capability as well as sterility). Such changes of material properties may in turn lead to treatment impairments, such as inadequate wound healing as well as infections. In the case of known allergies against components of SupraSDRM®, the membrane should not be applied. SupraSDRM® should be removed immediately if there are any signs of allergic reactions to the material. SupraSDRM® should be removed in cases of severe pain or accumulations of wound secretions. Coverage of intact skin may lead to skin macerations and should be avoided.


Additional Information

SupraSDRM® is gamma sterilised. SupraSDRM® is non-toxic. To be applied by healthcare professionals only. SupraSDRM® is a white membrane. In the case of a transparent membrane appearance please do not apply. As SupraSDRM® is a biosynthetic product, slight variations in shape and size may occur. These variations can be compensated by slight manual stretching of the membrane, if required. The summary of safety and clinical performance for SupraSDRM® is publicly available in the European Medical Device Database (Eudamed: <https://ec.europa.eu/tools/eudamed>) via the Basic UDI-DI 426018402AAA0000001PQ or via our homepage (<https://polymedics.com>).

Notice: any serious incident that has occurred in relation to SupraSDRM® should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. Users should inform the patients of residual risks associated with the medical device, particularly with respect to potentially serious infections, as indicated in these instructions for use at section contraindications and warnings and precautions.

This IFU applies to the REF numbers of SupraSDRM®, including 31D012-MDR, 31D018-MDR, 31D024-MDR, 310101-MDR, 310202-MDR, 310404-MDR, 310505-MDR, 310909-MDR, 310912-MDR, 311809-MDR, 311818-MDR.

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PolyMedics Innovations GmbH (PMI)