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| Caighdeáin Chomhréitithe | Cur Síos Gearr |
|------------------------------|--|
| Common specifications | Non available for the product |
| DIN EN ISO 13485 | Medical devices - Quality management systems - Requirements for regulatory purposes |
| DIN EN 62366-1 | Medical devices - Part 1: Application of usability engineering to medical devices |
| DIN EN ISO 14971 | Medical devices – Application of risk management to medical devices |
| DIN EN ISO 14155 | Clinical Investigation of Medical Devices for Human Subjects - Good Clinical Practice |
| DIN EN ISO 10993-1 | Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management system |
| DIN EN ISO 10993-3 | Biological evaluation of medical devices - Part 3: Tests for genotoxicity, carcinogenicity and reproductive toxicity |
| DIN EN ISO 10993-5 | Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity |
| DIN EN ISO 10993-6 | Biological evaluation of medical devices – Part 6: Tests for local effects after implantation |
| DIN EN ISO 10993-10 | Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization |
| DIN EN ISO 10993-11 | Biological evaluation of medical devices - Part 11: Tests for systemic toxicity |
| DIN EN ISO 10993-12 | Biological evaluation of medical devices - Part 12: Sample preparation and reference materials |
| DIN EN ISO 11737-1 | Sterilization of medical devices - Requirements for the estimation of population of microorganisms on a product |
| DIN EN ISO 11737-2 | Sterilization of medical devices - Microbiological methods - Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process |
| DIN EN ISO 11137-1 | Sterilization of health care products - Radiation - Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices |
| DIN EN ISO 11137-2 | Sterilization of Health Care Products - Radiation - Part 2: Establishing The Sterilization Dose |
| DIN EN 556-1 | Sterilization of medical devices - Requirements for medical devices to be designated "STERILE" - Part 1: Requirements for terminally sterilized medical devices |
| DIN EN ISO 11607-1 | Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems |

| | |
|---------------------------|--|
| DIN EN ISO 11607-2 | Packaging for terminally sterilized medical devices - Part 2: Validation requirements for forming, sealing and assembly processes |
| DIN EN ISO 20417 | Medical devices – Information to be supplied by the manufacturer |
| DIN EN ISO 15223-1 | Medical Devices - Symbols To Be Used With Medical Device Labels, Labelling And Information To Be Supplied - Part 1: General Requirements |
| DIN EN 868-2 | Packaging for terminally sterilized medical devices - Part 2: Sterilization wrap - Requirements and test methods |
| DIN EN 868-5 | Packaging for terminally sterilized medical devices - Part 5: Sealable pouches and reels of porous materials and plastic film construction - Requirements and test methods |
| DIN EN ISO 14698-1 | Cleanrooms and associated controlled environments -- Biocontamination control -- Part 1: General principles and methods |
| DIN EN ISO 14698-2 | Cleanrooms and associated controlled environments -- Biocontamination control -- Part 2: Evaluation and interpretation of biocontamination data |
| ISTA 2a | Partial Simulation Performance Tests - Packaged Products 150 lb (68 kg) or Less |
| USP <151> | Pyrogen Study |
| ASTM F1886/F1886M | Standard Test Method for Determining Integrity of Seals for Flexible Packaging by Visual Inspection |
| ASTM F88/F88M | Standard Test Method for Seal Strength of Flexible Barrier Materials |
| ASTM F3039 | Standard Test Method for Detecting Leaks in Nonporous Packaging or Flexible Barrier Materials by Dye Penetration |

9. Tagairtí sa litríocht

Baartmans, M. G., Dokter, J., Den Hollander, J. C., Kroon, A. A. & Oranje, A. P. (2011). Use of skin substitute dressings in the treatment of staphylococcal scalded skin syndrome in neonates and young infants. *Neonatology*, 100, 9-13.

Blome-Eberwein, S.A., Amani, H., Lozano, D.D., Gogal, C., Boorse, D., Pagella, P. Burns. A bio-degradable synthetic membrane to treat superficial and deep second degree burn wounds in adults and children – 4 year experience. 2020 Aug 29;S0305-4179(20)30507-6

Demircan, M., Gürünlüoğlu, K., Gözükar, Bag, H.G., Koçbıyık, A., Gül, M., Üremiş, N., Gül, S., Gürünlüoğlu, S., Türköz, Y., Taşçı, A. Impaction of the polylactic membrane or hydrofiber with silver dressings on the interleukin-6, tumor necrosis factor- α , transforming growth factor-b3 levels in the blood and tissues of pediatric patients with burns. *Ulus Travma Acil Cerrahi Derg.* 2021 Jan;27(1):122-131

Everett, M., Massand, S., Davis, W., Burkey, B. & Glat, P. M. (2015). Use of a copolymer dressing on superficial and partial-thickness burns in a paediatric population. *J Wound Care*, 24, S4-8.

Fischer, S., Kremer, T., Horter, J., Schaefer, A., Ziegler, B., Kneser, U., Hirche, C. Suprathel® for severe burns in the elderly: Case report and review of the literature. *Burns*. 2016 Aug;42(5):e86-92

Gürünlüoğlu, K., Demircan, M., Tasci, A., Uremis, M. M., Turkoz, Y., Bag, H. G., Akinci, A. & Bayrakci, E. (2019). The Effects of Two Different Burn Dressings on Serum Oxidative Stress Indicators in Children with Partial Burn. *J Burn Care Res*, 40, 444-450.

Gürünlüoğlu, K., Demircan, M., Koç, A., Koçbıyık, A., Taşçı, A., Durmuş, K., Gürünlüoğlu, S., Gözükara Bağ, H. The Effects of Different Burn Dressings on Length of Telomere and Expression of Telomerase in Children With Thermal Burns. *J Burn Care Res.* 2019 Apr 26;40(3):302-311.

Hakkarainen, T., Koivuniemi, R., Kosonen, M., Escobedo-Lucea, C., Sanz-Garcia, A., Vuola, J., Valtonen, J., Tammela, P., Mäkitie, A., Luukko, K., Yliperttula, M., Kavola, H. *J Control Release.* Nanofibrillar cellulose wound dressing in skin graft donor site treatment. 2016 Dec 28;244(Pt B):292-301.

Harenberg, P. S., Hrabowski, M., Ryssel, H., Gazyakan, E., Germann, G., Engel, H. & Reichenberger, M. A. (2010). CASE REPORT Febrile Ulceronecrotic Mucha-Habermann Disease. *Eplasty*, 10.

Highton, L., Wallace, C., Shah, M. Use of Suprathel® for partial thickness burns in children. *Burns.* 2013 Feb;39(1):136-41

Hundeshagen, G., Collins, V. N., Wurzer, P., Sherman, W., Voigt, C. D., Cambiaso-Daniel, J., Nunez Lopez, O., Sheaffer, J., Herndon, D. N., Finnerty, C. C. & Branski, L. K. (2018). A Prospective, Randomized, Controlled Trial Comparing the Outpatient Treatment of Pediatric and Adult Partial-Thickness Burns with Suprathel or Mepilex Ag. *J Burn Care Res*, 39, 261-267

Kaartinen, I. S. & Kuokkanen, H. O. (2011). Suprathel((R)) causes less bleeding and scarring than Mepilex((R)) Transfer in the treatment of donor sites of split-thickness skin grafts. *J Plast Surg Hand Surg*, 45, 200-3.

Keck, M., Selig, H. F., Lumenta, D. B., Kamolz, L. P., Mittlbock, M. & Frey, M. (2012). The use of Suprathel((R)) in deep dermal burns: first results of a prospective study. *Burns*, 38, 388-95.

Lindford, A. J., Kaartinen, I. S., Virolainen, S. & Vuola, J. (2011). Comparison of Suprathel(R) and allograft skin in the treatment of a severe case of toxic epidermal necrolysis. *Burns*, 37, e67-72.

Liodaki, E., Schopp, B.E., Lindert, J., Krämer, R., Kisch, T., Mailänder, P., Stang, F. Kombination von universellem Antidot und temporärem Hautersatz bei Verätzungen [Combination of a universal antidote and temporary skin substitute for chemical burns: Extended case report]. *Unfallchirurg.* 2015 Sep;118(9):804-7.

Madry, R., Struzyna, J., Stachura-Kulach, A., Drozd, L. & Bugaj, M. (2011). Effectiveness of Suprathel(R) application in partial thickness burns, frostbites and Lyell syndrome treatment. *Pol Przegl Chir*, 83, 541-8.

Markl, P., Prantl, L., Schreml, S., Babilas, P., Landthaler, M. & Schwarze, H. (2010). Management of split-thickness donor sites with synthetic wound dressings: results of a comparative clinical study. *Ann Plast Surg*, 65, 490-6.

Merz, K. M., Sievers, R., Reichert, B. (2011). Suprathel® for coverage of superficial dermal burns of the face. *GMS Verbrennungsmedizin*, 4.

Mueller, E., Haim, M., Petnehazy, T., Acham-Roschitz, B. & Trop, M. (2010). An innovative local treatment for staphylococcal scalded skin syndrome. *Eur J Clin Microbiol Infect Dis*, 29, 893-7.

Pfurtscheller, K., Zobel, G., Roedl, S. & Trop, M. (2008). Use of Suprathel dressing in a young infant with TEN. *Pediatr Dermatol*, 25, 541-3.

Pfurtscheller, K. & Trop, M. (2014). Phototoxic plant burns: report of a case and review of topical wound treatment in children. *Pediatr Dermatol*, 31, e156-9.

Rahmanian-Schwarz, A., Beiderwieden, A., Willkomm, L.M., Amr, A., Schaller, H.E., Lotter, O. A clinical evaluation of Biobrane[®] and Suprathel[®] in acute burns and reconstructive surgery. *Burns*. 2011 Dec;37(8):1343-8

Rashaan, Z. M., Krijnen, P., Allema, J. H., Vloemans, A. F., Schipper, I. B. & Breederveld, R. S. (2017). Usability and effectiveness of Suprathel[®] in partial thickness burns in children. *Eur J Trauma Emerg Surg*, 43, 549-556.

Rothenberger, J., Constantinescu, M. A., Held, M., Aebersold, D. M., Stolz, A., Tschumi, C. & Olariu, R. (2016). Use of a Polylactide-based Copolymer as a Temporary Skin Substitute for a Patient With Moist Desquamation Due to Radiation. *Wounds*, 28, E26-30.

Sari, E., Eryilmaz, T., Tetik, G., Ozakpinar, H. R. & Eker, E. (2014). Suprathel[®] -assisted surgical treatment of the hand in a dystrophic epidermolysis bullosa patient. *Int Wound J*, 11, 472-5.

Schiefer, J.L., Rahmanian-Schwarz, A., Schaller, H.E., Manoli, T. A Novel Hand-shaped Suprathel simplifies the Treatment of Partial-Thickness Burns. *Adv Skin Wound Care*. 2014 Nov;27(11):513-6

Schwarze, H., Kuntscher, M., Uhlig, C., Hierlemann, H., Prantl, L., Noack, N. & Hartmann, B. (2007). Suprathel, a new skin substitute, in the management of donor sites of split-thickness skin grafts: results of a clinical study. *Burns*, 33, 850-4.

Schwarze, H., Kuntscher, M., Uhlig, C., Hierlemann, H., Prantl, L., Ottomann, C. & Hartmann, B. (2008). Suprathel, a new skin substitute, in the management of partial-thickness burn wounds: results of a clinical study. *Ann Plast Surg*, 60, 181-5.

Selig, H. F., Keck, M., Lumenta, D. B., Mittlbock, M. & Kamolz, L. P. (2013). The use of a polylactide-based copolymer as a temporary skin substitute in deep dermal burns: 1-year follow-up results of a prospective clinical noninferiority trial. *Wound Repair Regen*, 21, 402-9.

Uhlig, C., Hierlemann, H., Dittel, K.-K. (2007). Actual Strategies in the Treatment of Severe Burns - Considering Modern Skin Substitutes. *Osteo trauma care* 15, 2-7.

Uhlig, C., Rapp, M. & Dittel, K. K. (2007a). [New strategies for the treatment of thermally injured hands with regard to the epithelial substitute Suprathel]. *Handchir Mikrochir Plast Chir*, 39, 314-9.

Uhlig, C., Rapp, M., Hartmann, B., Hierlemann, H., Planck, H. & Dittel, K. K. (2007b). Suprathel-an innovative, resorbable skin substitute for the treatment of burn victims. *Burns*, 33, 221-9.

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