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
Yhdenmukaistetut standardit	Lyhyt kuvaus
<b>Common specifications</b>	Non available for the product
<b>DIN EN ISO 13485</b>	Medical devices - Quality management systems - Requirements for regulatory purposes
<b>DIN EN 62366-1</b>	Medical devices - Part 1: Application of usability engineering to medical devices
<b>DIN EN ISO 14971</b>	Medical devices – Application of risk management to medical devices
<b>DIN EN ISO 14155</b>	Clinical Investigation of Medical Devices for Human Subjects - Good Clinical Practice
<b>DIN EN ISO 10993-1</b>	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management system
<b>DIN EN ISO 10993-3</b>	Biological evaluation of medical devices - Part 3: Tests for genotoxicity, carcinogenicity and reproductive toxicity

## Summary of Safety and Clinical Performance (SSCP)

<b>DIN EN ISO 10993-5</b>	Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity
<b>DIN EN ISO 10993-6</b>	Biological evaluation of medical devices – Part 6: Tests for local effects after implantation
<b>DIN EN ISO 10993-10</b>	Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization
<b>DIN EN ISO 10993-11</b>	Biological evaluation of medical devices - Part 11: Tests for systemic toxicity
<b>DIN EN ISO 10993-12</b>	Biological evaluation of medical devices - Part 12: Sample preparation and reference materials
<b>DIN EN ISO 11737-1</b>	Sterilization of medical devices - Requirements for the estimation of population of microorganisms on a product
<b>DIN EN ISO 11737-2</b>	Sterilization of medical devices - Microbiological methods - Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process
<b>DIN EN ISO 11137-1</b>	Sterilization of health care products - Radiation - Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices
<b>DIN EN ISO 11137-2</b>	Sterilization of Health Care Products - Radiation - Part 2: Establishing The Sterilization Dose
<b>DIN EN 556-1</b>	Sterilization of medical devices - Requirements for medical devices to be designated "STERILE" - Part 1: Requirements for terminally sterilized medical devices
<b>DIN EN ISO 11607-1</b>	Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems
<b>DIN EN ISO 11607-2</b>	Packaging for terminally sterilized medical devices - Part 2: Validation requirements for forming, sealing and assembly processes
<b>DIN EN ISO 20417</b>	Medical devices – Information to be supplied by the manufacturer
<b>DIN EN ISO 15223-1</b>	Medical Devices - Symbols To Be Used With Medical Device Labels, Labelling And Information To Be Supplied - Part 1: General Requirements
<b>DIN EN 868-2</b>	Packaging for terminally sterilized medical devices - Part 2: Sterilization wrap - Requirements and test methods
<b>DIN EN 868-5</b>	Packaging for terminally sterilized medical devices - Part 5: Sealable pouches and reels of porous materials and plastic film construction - Requirements and test methods
<b>DIN EN ISO 14698-1</b>	Cleanrooms and associated controlled environments -- Biocontamination control -- Part 1: General principles and methods
<b>DIN EN ISO 14698-2</b>	Cleanrooms and associated controlled environments -- Biocontamination control -- Part 2: Evaluation and interpretation of biocontamination data
<b>ISTA 2a</b>	Partial Simulation Performance Tests - Packaged Products 150 lb (68 kg) or Less
<b>USP &lt;151&gt;</b>	Pyrogen Study
<b>ASTM F1886/F1886M</b>	Standard Test Method for Determining Integrity of Seals for Flexible Packaging by Visual Inspection
<b>ASTM F88/F88M</b>	Standard Test Method for Seal Strength of Flexible Barrier Materials
<b>ASTM F3039</b>	Standard Test Method for Detecting Leaks in Nonporous Packaging or Flexible Barrier Materials by Dye Penetration

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Abbott, C. A., Carrington, A. L., Ashe, H., Bath, S., Every, L. C., Griffiths, J., Hann, A. W., Hussein, A., Jackson, N., Johnson, K. E., Ryder, C. H., Torkington, R., Van Ross, E. R., Whalley, A. M., Widdows, P.,

## Summary of Safety and Clinical Performance (SSCP)

Williamson, S., Boulton, A. J. & North-West Diabetes Foot Care, S. (2002). The North-West Diabetes Foot Care Study: incidence of, and risk factors for, new diabetic foot ulceration in a community-based patient cohort. *Diabet Med*, 19, 377-84.

Albertsson, A., Eklund, M. (1995). Influence of Molecular Structure on the degradation mechanism of Degradable Polymers: In Vitro Degradation of Poly(Trimethylene Carbonate), Poly(Trimethylene Carbonate-co-Caprolactone), and Poly(Adipic Anhydride). *J Applied Polymer Sci* 57, 87-103.

Ali, S. A., Zhong, S. P., Doherty, P. J. & Williams, D. F. (1993). Mechanisms of polymer degradation in implantable devices. I. Poly(caprolactone). *Biomaterials*, 14, 648-56.

Apelqvist, J., Bakker, K., Van Houtum, W. H., Nabuurs-Franssen, M. H. & Schaper, N. C. (2000). International consensus and practical guidelines on the management and the prevention of the diabetic foot. International Working Group on the Diabetic Foot. *Diabetes Metab Res Rev*, 16 Suppl 1, S84-92.

Baartmans, M. G., Dokter, J., Den Hollander, J. C., Kroon, A. A. & Oranje, A. P. (2011). Use of skin substitute dressings in the treatment of staphylococcal scalded skin syndrome in neonates and young infants. *Neonatology*, 100, 9-13.

Barbachowska A, Korzeniowski T, Surowiecka A, Tomaka P, Bugaj-Tobiasz M, Łączyk M, Górecka Z, Chrapusta A, Strużyna J. The Effectiveness of an Alloplastic Epidermal Substitute in the Treatment of Burn Wounds in Children: A Comparative Clinical Study of Skin Substitutes and Silver and Paraffin Gauze Dressings. *J Clin Med*. 2024 Nov 28;13(23):7238.

Barbachowska A, Tomaka P, Surowiecka A, Łączyk M, Górecka Z, Stepniewski A, Chrapusta A, Sadowy R, Strużyna J, Korzeniowski T (2025). The Influence of Blood Parameters on the Adhesion of an Epidermal Substitute in the Treatment of Burn Wounds in Children. *J Clin Med*. 2025 Jun 29;14(13):4614.

Blome-Eberwein, S.A., Amani, H., Lozano, D.D., Gogal, C., Boorse, D., Pagella, P. Burns. A bio-degradable synthetic membrane to treat superficial and deep second degree burn wounds in adults and children – 4 year experience. 2020 Aug 29;S0305-4179(20)30507-6

Bostman, O. M. (1991). Absorbable implants for the fixation of fractures. *J Bone Joint Surg Am*, 73, 148-53.

Brady, J. M., Cutright, D. E., Miller, R. A. & Barristone, G. C. (1973). Resorption rate, route, route of elimination, and ultrastructure of the implant site of polylactic acid in the abdominal wall of the rat. *J Biomed Mater Res*, 7, 155-66.

Burd, A. & Yuen, C. (2005). A global study of hospitalized paediatric burn patients. *Burns*, 31, 432-8.

Bryant, R., Nix, D. (2006). *Acute and Chronic Wounds* 3rd Edition.

Cha, Y. & Pitt, C. G. (1990). The biodegradability of polyester blends. *Biomaterials*, 11, 108-12.

Cheema L, Manzoor S, Khalid U, Shamim R, Hashaam, Tayyab Z, Bashir M. Suprathel Dressing at Split Thickness Skin Graft Donor Site for Pain Control and Wound Healing. *Pakistan journal of medical and health sciences*, 2022, 16(10), 116-118

Dadras M, Wagner JM, Wallner C, Sogorski A, Sacher M, Harati K, Lehnhardt M, Behr B. Enzymatic debridement of hands with deep burns: a single center experience in the treatment of 52 hands. *J Plast Surg Hand Surg*. 2020 Aug;54(4):220-224. doi: 10.1080/2000656X.2020.1755298. Epub 2020 Apr 21. PMID: 32315252.

## Summary of Safety and Clinical Performance (SSCP)

Dastagir N, Kijas D, Obed D, Tamulevicius M, Vogt PM, Dastagir K. Suprathel® and water-filtered infrared-A radiation (wIRA) as a new treatment strategy for toxic epidermal necrolysis (TEN): A prospective study. *Burns*. 2024 Dec;50(9):107283

Delgado-Miguel C, García Morán A, Fuentes Gómez L, Díaz M, Miguel-Ferrero M, López-Gutiérrez JC. Comparison of the effectiveness of three different skin substitutes for the treatment of pediatric burns. *Eur J Pediatr*. 2024 Dec 13;184(1):80.

Demircan, M., Gürünlüoğlu, K., Gözü kara, Bag, H.G., Koçbıyık, A., Gül, M., Üremiş, N., Gül, S., Gürünlüoğlu, S., Türköz, Y., Taşçı, A. Impaction of the polylactic membrane or hydrofiber with silver dressings on the interleukin-6, tumor necrosis factor- $\alpha$ , transforming growth factor-b3 levels in the blood and tissues of pediatric patients with burns. *Ulus Travma Acil Cerrahi Derg*. 2021 Jan;27(1):122-131

Dhivya, S., Padma, V. V. & Santhini, E. (2015). Wound dressings - a review. *Biomedicine (Taipei)*, 5, 22.

Everett, M., Massand, S., Davis, W., Burkey, B. & Glat, P. M. (2015). Use of a copolymer dressing on superficial and partial-thickness burns in a paediatric population. *J Wound Care*, 24, S4-8.

Fernandes S, Teixeira I, Carmo L, Campos M, Garcia M. The Use of a Polylactic Membrane in Pediatric Burns Proves to be Successful Even After Late Application. *J Burn Care Res*. . 2023 Sep 7;44(5):1176-1181

Fischer, S., Kremer, T., Horter, J., Schaefer, A., Ziegler, B., Kneser, U., Hirche, C. Suprathel® for severe burns in the elderly: Case report and review of the literature. *Burns*. 2016 Aug;42(5):e86-92

Galati V, Vonthein R, Stang F, Mailaender P, Kisch T. Split thickness skin graft versus application of the temporary skin substitute suprathel in the treatment of deep dermal hand burns: a retrospective cohort study of scar elasticity and perfusion. *Int J Burns Trauma*. 2021 Aug 15;11(4):312-320

Gürünlüoğlu, K., Demircan, M., Tasci, A., Uremis, M. M., Turkoz, Y., Bag, H. G., Akinci, A. & Bayrakci, E.. The Effects of Two Different Burn Dressings on Serum Oxidative Stress Indicators in Children with Partial Burn. *J Burn Care Res*, 40, 444-450 (2019a).


Gürünlüoğlu, K., Demircan, M., Koç, A., Koçbıyık, A., Taşçı, A., Durmuş, K., Gürünlüoğlu, S., Gözü kara Bağ, H. The Effects of Different Burn Dressings on Length of Telomere and Expression of Telomerase in Children With Thermal Burns. *J Burn Care Res*. 2019 Apr 26;40(3):302-311 (2019b).

Hakkarainen, T., Koivuniemi, R., Kosonen, M., Escobedo-Lucea, C., Sanz-Garcia, A., Vuola, J., Valtonen, J., Tammela, P., Mäkitie, A., Luukko, K., Yliperttula, M., Kavola, H. J Control Release. Nanofibrillar cellulose wound dressing in skin graft donor site treatment. 2016 Dec 28;244(Pt B):292-301.

Harenberg, P. S., Hrabowski, M., Rysse, H., Gazyakan, E., Germann, G., Engel, H. & Reichenberger, M. A. (2010). CASE REPORT Febrile Ulceronecrotic Mucha-Habermann Disease. *Eplasty*, 10.

Heitzmann W, Enzmann J, von Kohout M, Mattern MM, Akkan J, Fuchs PC, Schiefer JL. Accelerated wound healing of enzymatically debrided deep dermal burn wounds after the use of fish skin (Kerecis Omega3 Wound®) in comparison to Suprathel®. *Burns*. 2025 Jun;51(5):107471. (2025a)

Heitzmann W, Mossing M, Fuchs PC, Akkan J, Seyhan H, Grieb G, Opländer C, Schiefer JL. Comparative Clinical Study of Suprathel® and Jelonet® Wound Dressings in Burn Wound Healing after Enzymatic Debridement.

## Summary of Safety and Clinical Performance (SSCP)

Heitzmann W, Mossing M, von Kohout M, Alhujaili Y, Akkan J, Fuchs PC, Lefering R, Schiefer JL. Long-term skin quality and scar formation after enzymatic debridement of deep-dermal burn wounds: A follow-up comparative study of Suprathel® and Jelonet®. *Burns*. 2025 Nov 20;52(1):107803. (2025b)

Hettiaratchy, S. & Papini, R. (2004). Initial management of a major burn: II--assessment and resuscitation. *BMJ*, 329, 101-3.

Highton, L., Wallace, C., Shah, M. Use of Suprathel® for partial thickness burns in children. *Burns*. 2013 Feb;39(1):136-41

Hollinger, J. O. & Battistone, G. C. (1986). Biodegradable bone repair materials. Synthetic polymers and ceramics. *Clin Orthop Relat Res*, 290-305.

Hundeshagen, G., Collins, V. N., Wurzer, P., Sherman, W., Voigt, C. D., Cambiaso-Daniel, J., Nunez Lopez, O., Sheaffer, J., Herndon, D. N., Finnerty, C. C. & Branski, L. K. (2018). A Prospective, Randomized, Controlled Trial Comparing the Outpatient Treatment of Pediatric and Adult Partial-Thickness Burns with Suprathel or Mepilex Ag. *J Burn Care Res*, 39, 261-267.

Hunt, T., Aslam, R., Beckert, S., Wagner, S., Ghani, R., Hussain, M., Roy, S., And Sen, C. (2007). Aerobically Derived Lactate Stimulates Revascularization And Tissue Repair Via Redox Mechanisms. *Antioxidants & Redox Signaling* Volume 9, Number 8.

Kaartinen, I. S. & Kuokkanen, H. O. (2011). Suprathel((R)) causes less bleeding and scarring than Mepilex((R)) Transfer in the treatment of donor sites of split-thickness skin grafts. *J Plast Surg Hand Surg*, 45, 200-3.

Karlsson M, Steinvall I, Elmasry M. Suprathel® or Mepilex® Ag for treatment of partial thickness burns in children: A case control study. *Burns*. 2023 Mar 11;S0305-4179(23)00043-8.

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Kamolz, L., Lumenta, D., Kitzinger, H., Frey, M. (2008). Tissue engineering for cutaneous wounds: an overview of current standards and possibilities. *Eur Surg* 40, 19-26.

Katz, A. R., Mukherjee, D. P., Kaganov, A. L. & Gordon, S. (1985). A new synthetic monofilament absorbable suture made from polytrimethylene carbonate. *Surg Gynecol Obstet*, 161, 213-22.

Keck, M., Selig, H. F., Lumenta, D. B., Kamolz, L. P., Mittlbock, M. & Frey, M. (2012). The use of Suprathel((R)) in deep dermal burns: first results of a prospective study. *Burns*, 38, 388-95.

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Lipkin I, Hughes M, Hughes WB (2025) Utilization of Suprathel® in delayed surgical management of toxic epidermal necrolysis (TEN): A case report, *Burns Open* 2025

Madry, R., Struzyna, J., Stachura-Kulach, A., Drozd, L. & Bugaj, M. (2011). Effectiveness of Suprathel(R) application in partial thickness burns, frostbites and Lyell syndrome treatment. *Pol Przegl Chir*, 83, 541-8.

Margolis, D. J., Allen-Taylor, L., Hoffstad, O. & Berlin, J. A. (2004). The accuracy of venous leg ulcer prognostic models in a wound care system. *Wound Repair Regen*, 12, 163-8.

Margolis, D. J., Allen-Taylor, L., Hoffstad, O. & Berlin, J. A. (2004). The accuracy of venous leg ulcer prognostic models in a wound care system. *Wound Repair Regen*, 12, 163-8.

Margolis, D. J., Malay, D. S., Hoffstad, O. J., Leonard, C. E., Macurdy, T., Lopez De Nava, K., Tan, Y., Molina, T. & Siegel, K. L. 2011. Prevalence of diabetes, diabetic foot ulcer, and lower extremity amputation among Medicare beneficiaries, 2006 to 2008: Data Points #1. *Data Points Publication Series*. Rockville (MD).

Markl, P., Prantl, L., Schreml, S., Babilas, P., Landthaler, M. & Schwarze, H. (2010). Management of split-thickness donor sites with synthetic wound dressings: results of a comparative clinical study. *Ann Plast Surg*, 65, 490-6.

Maxey J, Wampler M, Elkady D, DeVault A, Williamson K, Borrelli M, Lou R, Khandelwal A. Minimally invasive excision combined with epidermal autografting and poly-lactic acid skin substitute improves outcomes in pediatric partial thickness burns, *Journal of Pediatric Surgery Open*, Volume 11 (2025)

März V, Vogt M. Skin Healing of Deep Second Degree Burn Injuries in Four Individuals Sustained in a Boat Explosion Results after Different Approaches. *Eur. Burn J*. 2020, 1, 191–195


Merz, K. M., Sievers, R., Reichert, B. (2011). Suprathel® for coverage of superficial dermal burns of the face. *GMS Verbrennungsmedizin*, 4.

Miguel-Ferrero M, Delgado-Miguel C, Díaz M, Carlos López-Gutiérrez J. Toxic epidermal necrolysis management with suprathel™ Tratamiento de la necrólisis epidérmica tóxica con suprathel®. *An Pediatr (Engl Ed)*. 2023 Aug 17:S2341-2879(23)00185-0

Moellhoff N, Lettner M, Frank K, Giunta RE, Ehrl D. Poly-lactic Acid Membrane Improves Outcome of Split-Thickness Skin Graft Donor Sites: A Prospective, Comparative, Randomized Study. *Plast Reconstr Surg*. 2022 Nov 1;150(5):1104-1113

Moffatt, C. J. & Dorman, M. C. (1995). Recurrence of leg ulcers within a community ulcer service. *J Wound Care*, 4, 57-61.

Moher D, Liberati A, Tetzlaff J, Altman DG, The PRISMA Group. Preferred Reporting Items for Systematic Reviews and Meta-Analyses: The PRISMA Statement. *PLoS Med*. 2009;6:e1000097. DOI:10.1371/journal.pmed1000097

## Summary of Safety and Clinical Performance (SSCP)

Mueller, E., Haim, M., Petnehazy, T., Acham-Roschitz, B. & Trop, M. (2010). An innovative local treatment for staphylococcal scalded skin syndrome. *Eur J Clin Microbiol Infect Dis*, 29, 893-7.

Novoa-Moreno AL, Martinez-Jimenez MA, Ortiz-Alvarez A, Sanchez-Olivo N, Loza-Gonzalez VM, Ramirez-GarciaLuna JL (2025). Polylactic Acid Membranes, a Novel Adjunct Treatment for Bullous Impetigo. *Infect Dis Rep*. 2025 Jun 19;17(3):72.

Nischwitz SP, Popp D, Shubitidze D, Luze H, Zrim R, Klemm K, Rapp M, Haller HL, Feisst M, Kamolz LP. The successful use of polylactide wound dressings for chronic lower leg wounds: A retrospective analysis. *Int Wound J*. 2021 Nov 8.

O'meara, S. & Martyn-St James, M. (2013). Foam dressings for venous leg ulcers. *Cochrane Database Syst Rev*, CD009907

Pego, A. P., Van Luyn, M. J., Brouwer, L. A., Van Wachem, P. B., Poot, A. A., Grijpma, D. W. & Feijen, J. (2003). In vivo behavior of poly(1,3-trimethylene carbonate) and copolymers of 1,3-trimethylene carbonate with D,L-lactide or epsilon-caprolactone: Degradation and tissue response. *J Biomed Mater Res A*, 67, 1044-54.

Pfurtscheller, K. & Trop, M. (2014). Phototoxic plant burns: report of a case and review of topical wound treatment in children. *Pediatr Dermatol*, 31, e156-9.

Pfurtscheller, K., Zobel, G., Roedl, S. & Trop, M. (2008). Use of Suprathel dressing in a young infant with TEN. *Pediatr Dermatol*, 25, 541-3.

Pitt, C. G., Gratzl, M. M., Kimmel, G. L., Surles, J. & Schindler, A. (1981). Aliphatic polyesters II. The degradation of poly (DL-lactide), poly (epsilon-caprolactone), and their copolymers in vivo. *Biomaterials*, 2, 215-20.

Quinn, K. J., Courtney, J. M., Evans, J. H., Gaylor, J. D. & Reid, W. H. (1985). Principles of burn dressings. *Biomaterials*, 6, 369-77.

Rahmanian-Schwarz, A., Beiderwieden, A., Willkomm, L.M., Amr, A., Schaller, H.E., Lotter, O. A clinical evaluation of Biobrane® and Suprathel® in acute burns and reconstructive surgery. *Burns*. 2011 Dec;37(8):1343-8

Rajendran, S., Anand, S.C. (2011). Hi-tech textiles for interactive wound therapies: Handbook of Medical Textiles.


Rashaan, Z. M., Krijnen, P., Allema, J. H., Vloemans, A. F., Schipper, I. B. & Breederveld, R. S. (2017). Usability and effectiveness of Suprathel® in partial thickness burns in children. *Eur J Trauma Emerg Surg*, 43, 549-556.

Raine BE, Fowler CC, Nguyen A, Schneider P, Gunn KJ, Chikoti R, Bell DE (2025). Polylactic Acid Membrane for the Treatment of Pediatric Burn Injuries: An Adoptable Practice to Improve Burn Management. *Ann Plast Surg*. 2025 Apr 1;94(4S Suppl 2)

Robson, M. C., Steed, D. L. & Franz, M. G. (2001). Wound healing: biologic features and approaches to maximize healing trajectories. *Curr Probl Surg*, 38, 72-140.

Rothenberger, J., Constantinescu, M. A., Held, M., Aebersold, D. M., Stolz, A., Tschumi, C. & Olariu, R. (2016). Use of a Polylactide-based Copolymer as a Temporary Skin Substitute for a Patient With Moist Desquamation Due to Radiation. *Wounds*, 28, E26-30.

Ruckley, C. V. (1998). Caring for patients with chronic leg ulcer. *BMJ*, 316, 407-8.

## Summary of Safety and Clinical Performance (SSCP)

- Sari, E., Eryilmaz, T., Tetik, G., Ozakpinar, H. R. & Eker, E. (2014). Suprathel((R)) -assisted surgical treatment of the hand in a dystrophic epidermolysis bullosa patient. *Int Wound J*, 11, 472-5.
- Schiefer, J.L., Rahmanian-Schwarz, A., Schaller, H.E., Manoli, T. A Novel Hand-shaped Suprathel simplifies the Treatment of Partial-Thickness Burns. *Adv Skin Wound Care*. 2014 Nov;27(11):513-6
- Schiefer JL, Andreae J, Bagheri M, Fuchs PC, Lefering R, Heitzmann W, Schulz A (2022a). A clinical comparison of pure knitted silk and a complex synthetic skin substitute for the treatment of partial thickness burns. *Int Wound J*. 2022 Jan;19(1):178-187.
- Schiefer JL, Aretz GF, Fuchs PC, Bagheri M, Funk M, Schulz A, Daniels M (2022b). Comparison of wound healing and patient comfort in partial-thickness burn wounds treated with SUPRATHEL and epicte hydro wound dressings. *Int Wound J*. 2022 May;19(4):782-790.
- Schiefer JL, Aretz FG, Fuchs PC, Lefering R, Yary P, Opländer C, Schulz A, Daniels M (2022c). Comparison of Long-Term Skin Quality and Scar Formation in Partial-Thickness Burn Wounds Treated with Suprathel® and epicitehydro® Wound Dressings. *Medicina (Kaunas)*. 2022 Oct 28;58(11):1550.
- Schiefer JL, Andreae J, Fuchs PC, Lefering R, Heidekrueger PI, Schulz A, Bagheri M (2022d). Evaluation of Scar Quality after Treatment of Superficial Burns with Dressilk® and Suprathel®-In an Intraindividual Clinical Setting. *J Clin Med*. 2022 May 18;11(10):2857
- Schreml, S., Szeimies, R. M., Prantl, L., Karrer, S., Landthaler, M. & Babilas, P. (2010). Oxygen in acute and chronic wound healing. *Br J Dermatol*, 163, 257-68.
- Schriek K, Ott H, Sinnig M. Paradigm Shift in Treatment Strategies for Second-Degree. *Eur. Burn J*. 2022, 3, 1-9
- Schulz A, Perbix W, Shoham Y, Daali S, Charalampaki C, Fuchs PC, Schiefer J. Our initial learning curve in the enzymatic debridement of severely burned hands-Management and pit falls of initial treatments and our development of a post debridement wound treatment algorithm. *Burns*. 2017 Mar;43(2):326-336. doi: 10.1016/j.burns.2016.08.009. Epub 2016 Oct 27. PMID: 28341257.
- Schwarze, H., Kuntscher, M., Uhlig, C., Hierlemann, H., Prantl, L., Noack, N. & Hartmann, B. (2007). Suprathel, a new skin substitute, in the management of donor sites of split-thickness skin grafts: results of a clinical study. *Burns*, 33, 850-4.
- Schwarze, H., Kuntscher, M., Uhlig, C., Hierlemann, H., Prantl, L., Ottomann, C. & Hartmann, B. (2008). Suprathel, a new skin substitute, in the management of partial-thickness burn wounds: results of a clinical study. *Ann Plast Surg*, 60, 181-5.
- Selig, H. F., Keck, M., Lumenta, D. B., Mittlbock, M. & Kamolz, L. P. (2013). The use of a polylactide-based copolymer as a temporary skin substitute in deep dermal burns: 1-year follow-up results of a prospective clinical noninferiority trial. *Wound Repair Regen*, 21, 402-9.
- Sun, H., Mei, L., Song, C., Cui, X. & Wang, P. (2006). The in vivo degradation, absorption and excretion of PCL-based implant. *Biomaterials*, 27, 1735-40.
- Szycher, M. & Lee, S. J. (1992). Modern wound dressings: a systematic approach to wound healing. *J Biomater Appl*, 7, 142-213.
- Tams, J., Joziassse, C. A., Bos, R. R., Rozema, F. R., Grijpma, D. W. & Pennings, A. J. (1995). High-impact poly(L/D-lactide) for fracture fixation: in vitro degradation and animal pilot study. *Biomaterials*, 16, 1409-15.

## Summary of Safety and Clinical Performance (SSCP)

Thomas, S. S., Lawrence, J. C. & Thomas, A. (1995). Evaluation of hydrocolloids and topical Medication In Minor Burns. *J Wound Care*, 4, 218-20.

Trabold, O., Wagner, S., Wicke, C., Scheuenstuhl, H., Hussain, Z., Rosen, N., Seremetiev, A., Becker, H., Hunt, T (2003). Lactate and oxygen constitute a fundamental regulatory mechanism in wound healing. *Wound Rep Reg* 11:504–509.

Uhlig, C., Hierlemann, H., Dittel, K.-K. (2007). Actual Strategies in the Treatment of Severe Burns - Considering Modern Skin Substitutes. *Osteo trauma care* 15, 2-7.

Uhlig, C., Rapp, M. & Dittel, K. K. (2007a). [New strategies for the treatment of thermally injured hands with regard to the epithelial substitute Suprathel]. *Handchir Mikrochir Plast Chir*, 39, 314-9.

Uhlig, C., Rapp, M., Hartmann, B., Hierlemann, H., Planck, H. & Dittel, K. K. (2007b). Suprathel-an innovative, resorbable skin substitute for the treatment of burn victims. *Burns*, 33, 221-9.

van de Warenburg MS, El Yadari S, Hummelink S, Ulrich DJ, Vehmeijer-Heeman M (2025b). Suprathel's usability and effectiveness for the treatment of paediatric partial thickness burns: a 10-year retrospective cohort study. *Burns*. 2025 May;51(4):107451.

van de Warenburg MS, Teeuwen B, Hummelink S, Ulrich DJ, Vehmeijer-Heeman ML (2025a). Does the dressing matter in pediatric partial-thickness burns: a systematic review and meta-analysis. *Burns*. 2025 May;51(4):107428.

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Vowden, K. R. & Vowden, P. (2009). The prevalence, management, equipment provision and outcome for patients with pressure ulceration identified in a wound care survey within one English health care district. *J Tissue Viability*, 18, 20-6.

Wallner C, Holtermann J, Drysch M, Schmidt S, Reinkemeier F, Wagner JM, Dadras M, Sogorski A, Houschyar KS, Becerikli M, Lehnhardt M, Behr B. The Use of Intact Fish Skin as a Novel Treatment Method for Deep Dermal Burns Following Enzymatic Debridement: A Retrospective Case-Control Study. *Eur Burn J*. 2022 Jan 27;3(1):43-55. doi: 10.3390/ebj3010006. PMID: 39604176; PMCID: PMC11575382.

Wasiak, J., Cleland, H., Campbell, F. & Spinks, A. (2013). Dressings for superficial and partial thickness burns. *Cochrane Database Syst Rev*, CD002106.

Wasiak, J., Cleland, H. (2005). Minor thermal burns. *Clinical Evidence*, 2754-2763.

Wild, S., Roglic, G., Green, A., Sicree, R. & King, H. (2004). Global prevalence of diabetes: estimates for the year 2000 and projections for 2030. *Diabetes Care*, 27, 1047-53.

Woodward, S. C., Brewer, P. S., Moatamed, F., Schindler, A. & Pitt, C. G. (1985). The intracellular degradation of poly(epsilon-caprolactone). *J Biomed Mater Res*, 19, 437-44.

Wu, L., Norman, G., Dumville, J. C., O'meara, S. & Bell-Syer, S. E. (2015). Dressings for treating foot ulcers in people with diabetes: an overview of systematic reviews. *Cochrane Database Syst Rev*, CD010471.

Yilmaz E, Simsek E, Ertenlice A, Akgun AE, Akin M (2025). Use of Suprathel for Deep Dermal Burns: Our Clinical Experience. *ACH Med J* (2025) 3: 089-092

## Summary of Safety and Clinical Performance (SSCP)

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