INSTRUCTIONS FOR USE

EN	Alloplastic skin substitute
DE	Alloplastischer Hautersatz
BG	Алопластичен кожен заместител
CS	Aloplastická kožní náhrada
DA	Alloplastisk huderstatning
EL	Αλλοπλαστικό υποκατάστατο δέρματος
ES	Reemplazo aloplástico de la piel
ET	Alloplastikas kasutatav naha aseaine
FI	Alloplastinen ihon korvike
FR	Substitut alloplastique de la peau
GA	Leaschraiceann allaplaisteach
HR	Aloplastični nadomjestak za kožu
HU	Alloplasztikus bőrhelyettesítő
IT	Sostituto alloplastico della cute
LT	Aloplastinis odos pakaitalas
LV	Aloplastisks ādas aizvietotājs
MT	Sostitut tal-ġilda alloplastiku
NL	Alloplastische huidvervanger
NO	Alloplastisk huderstatning
PL	Alloplastyczny substytut skóry
PT	Substituto aloplástico da pele
RO	Inlocuitor aloplastic al pielii
SK	Aloplastická náhrada kože
SL	Aloplastični nadomestek kože
SV	Alloplastiskt hudsubstitut
TR	Alloplastik deri ikamesi

EN Instructions for use

Product description

Suprathel® is an absorbable, microporous membrane and approved as alloplastic skin substitute for the treatment of epidermal and dermal wounds.

Composition

Copolymer composed out of lactide (73–80 %), trimethylene carbonate (12–16 %) and caprolactone (8–11 %).

Indications

Suprathel® is used for patients with epidermal and dermal wounds, such as abrasions, split skin graft donor sites, 2nd degree burns as well as 2nd degree burns mixed with 3rd degree burned areas. Suprathel® is used for patients with chronic wounds, such as venous and arterial ulcers, as well as diabetic wounds.

Contraindications

Suprathel® should not be used on infected wound sites or on severe bleeding wounds without additional hemostatic treatment. Suprathel® should not be applied on chronic dry wounds.

Properties and Mechanism

Suprathel® is an alloplastic, absorbable skin substitute that is highly permeable to oxygen and water vapour. Suprathel® can be adapted to any body part due to its elastic/plastic properties.

Application Procedure

Clean and debride the wound prior to applying Suprathel®. Remove Suprathel® from the package in an aseptic manner. Be sure to separate and discard the green colored protective paper that surrounds the white Suprathel® membrane. Apply Suprathel® on the affected area in either direction or side. If multiple dressings are required to cover larger areas, Suprathel® may be applied overlapping to assure complete coverage of the wound area. Slight stretching of Suprathel® may facilitate contouring of the dressing around curved anatomies such as elbows, fingers, etc. Assure complete contact of Suprathel® with the affected area by removing all air bubbles and wrinkles. Suprathel® may be trimmed so that only injured tissues are covered by the dressing. For protection against mechanical disruptions, apply a single sheet of fatty

gauze or similar dressing material over Suprathel®.

Aftercare Procedure

A gauze dressing applied over the fatty gauze can act as a protective barrier against mechanical alteration or contamination. This dressing can be changed as needed. Removal of Suprathel® during such dressing changes should be avoided to prevent disruption of the underlying healing tissues below the dressing. In case of chronic wounds, depending on the condition of the wound, it may be necessary to change or reapply Suprathel® at certain intervals. The transparent nature of applied Suprathel® facilitates visualization of healing tissue.

Adverse reactions

There have been no serious adverse side effects attributed to the product.

Storage and Handling

Suprathel® may only be stored in original packaging dry and cool (-20 up to +22 °C). Keep the product away from heat sources. Suprathel® can be disposed of with household waste.

Warnings and Precautions

Do not apply a product, where the sterility may not be ensured as this may lead to severe infections. The content is sterile unless sterile packaging is damaged. In case of packaging damages, the sterility of the product is not ensured. The unused contents of opened or damaged sterile packages are to be discarded.

Do not reuse and do not resterilise. If the product is nevertheless reused, this may lead to impairment of product performance characteristics (reduced permeability, elasticity, adherence capability as well as sterility).

Such changes of material properties may in turn lead to treatment impairments, such as inadequate wound healing as well as infections. In the case of known allergies against components of Suprathel®, the membrane should not be applied. Suprathel® should be removed immediately if there are any signs of allergic reactions to the material. Suprathel® should be removed in cases of severe pain or accumulations of wound secretions. Coverage of intact skin may lead to skin macerations and

should be avoided.

Additional Information


Suprathel® is gamma sterilised. Suprathel® is non-toxic. To be applied by healthcare professionals only. Suprathel® is a white membrane. In the case of a transparent membrane appearance please do not apply. As Suprathel® is a biosynthetic product, slight variations in shape and size may occur. These variations can be compensated by slight manual stretching of the membrane, if required.

The summary of safety and clinical performance for Suprathel® is publicly available in the European Medical Device Database (Eudamed: <https://ec.europa.eu/tools/eudamed>) via the Basic UDI-DI 426018402AAA0000001PQ or via our homepage (<https://polymedics.com>).

Notice: any serious incident that has occurred in relation to Suprathel® should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. Users should inform the patients of residual risks associated with the medical device, particularly with respect to potentially serious infections, as indicated in these instructions for use at section contraindications and warnings and precautions.

This IFU applies to the REF numbers of Suprathel®, including 110505-MDR, 150505-MDR, 110910-MDR, 150910-MDR, 111810-MDR, 151810-MDR, 111823-MDR, 151823-MDR, 11F001-MDR and 12H001-MDR.

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PolyMedics Innovations GmbH (PMI)