


**SUPRATHEL**<sup>®</sup>

**SKIN**

**THE TEMPORARY SECOND SKIN**



The Wound Healing Company



## Intelligent wound care with a temporary skin substitute

**SUPRATHEL®** is an innovative skin substitute indicated for the treatment of epidermal and dermal wounds. Successful use of the product has been demonstrated in the management of burns and STSG donor sites.

Just like a second skin, SUPRATHEL® covers the wound and leads the wound through a quick, complication free healing process. SUPRATHEL® was developed analogous to the human skin and thus shares the same properties such as elasticity, permeability to water vapor and impermeability to bacteria.

SUPRATHEL® is a single application product that is applied directly to a disinfected and debrided wound bed, where it stays intact until the wound is completely healed. After it is applied, the membrane becomes translucent and makes inspection of the healing process possible.

## Indications for use


**SUPRATHEL®** was developed for the treatment of epidermal and dermal wounds.

**SUPRATHEL®** covers a wide range of wound care:

- Split-thickness skin graft (STSG) donor sites
  
- Burns
  - Superficial
  - Partial thickness
  - Partial thickness with 3° areas



## Positioning in the treatment of wounds



\*Full thickness burns not indicated in the US and Canada


- Alginate
- Hydrofibres
- Hydrogels
- Foam dressing
- Hydrocolloids
- Film dressing
- Collagen dressing

- Cadaver-based scaffolds
- STSG
- Mesh-graft transplantations
- Cultured epithelial autografts (CEA)
- Acellular grafts
- Dermal substitutes
- Xenograft

## Properties

<b>Composition</b>	Lacto-capromer
<b>Degradation</b>	Hydrolytically
<b>Plasticity</b>	> 50 % elongation at break
<b>Permeability to water vapor</b>	40 - 70 ml/m <sup>2</sup> (hour) approx. 1.000 - 1.700 per day
<b>Porosity</b>	70 - 80 %

## The unique combination of strong features




\* The individual components have been successfully used in surgery for decades.

**SUPRATHEL® is purely synthetic and therefore does not bear any residual risks as is the case with biological products of human or animal origin.**

Literature suggests that lactate may stimulate the healing process by supporting angiogenesis<sup>1-6</sup> and the re-building of the dermis<sup>7-10</sup>. The potential of lactate to act as a free radical scavenger and therefore to be able to reduce oxidative stress has been demonstrated within the literature<sup>11</sup>.

SUPRATHEL® TEMPORARY SKIN SUBSTITUTE

DAMAGED SKIN PART



### Literature

- <sup>1</sup> Lu et al. 2002: J. Biol. Chem. 277:23111-5.
- <sup>2</sup> Lu et al. 2005: J. Biol. Chem. 280:41928-39.
- <sup>3</sup> Constant et al. 2000: Wound. Repair Regen. 8:353-360.
- <sup>4</sup> Rendl et al. 2001: Br. J. Dermatol. 145:3-9.
- <sup>5</sup> Beckert et al. 2006: Wound. Repair Regen. 14: 321-324.
- <sup>6</sup> Nareike et al. 2005: Am. J. Physiol. Endocrinol. Metab. 289:E534-42.

- <sup>7</sup> Green and Goldberg 1964: Nature 204: 347-9.
- <sup>8</sup> Hunt et al. 1978: Am J Surg. 135(3):328-32.
- <sup>9</sup> Klein et al. 2001: J Hand Surg Am 26(5):847-54.
- <sup>10</sup> Wagner et al. 2004: Wound. Repair Regen. 12:368-73.
- <sup>11</sup> Groussard et al. 2000: J Appl Physiol. 89: 169-175.



## The advantages of working with SUPRATHEL®

- **One-time wound dressing, no change of SUPRATHEL® needed**

- **Significant pain relief<sup>1-4</sup> – by up to 60%**

- » Less anesthesia and pain medication required<sup>1,5,6</sup>

- **Low rate of infections<sup>1,2,5,7-9</sup>, no biological risk**

- » Synthetic, biocompatible, absorbable
- » No reported allergic reactions

- **Fast wound healing<sup>1,2,5</sup>**

- » Improved epithelization<sup>3,11</sup>

- **Lower treatment costs<sup>2,4,5</sup>**

- » Less cost and effort for dressing changes<sup>2,9</sup>
- » Shorter hospital stay<sup>5</sup>

- **Excellent cosmetic outcomes and scar quality<sup>4,7,8</sup>**

- **Low inflammatory reaction<sup>12</sup>**

- **Reduced transplantation rate<sup>6</sup>**

### Literature

<sup>1</sup> Schwarze et al. 2007: Burns. 2007 Nov;33(7):850-4

<sup>2</sup> Schwarze et al. 2008: Ann Plast Surg. 2008 Feb;60(2):181-5

<sup>3</sup> Uhlig et al. 2007: Burns. 2007 Mar;33(2):221-9

<sup>4</sup> Kaartinen and Kuokkanen 2011: Burns. 2012 May;38(3):388-95.

<sup>5</sup> Everett et al. 2015: J Wound Care. 2015 Jul;24(7):S4-8.

<sup>6</sup> Schriek et al. 2022: Eur. Burn J. 2022, 3, 1–9

<sup>7</sup> Keck et al. 2012: Burns. 2012 May;38(3):388-95.

<sup>8</sup> Hundeshagen et al. 2018: J Burn Care Res. 2018 Feb 20;39(2):261-267.

<sup>9</sup> Markl et al. 2010: Ann Plast Surg. 2010 Nov;65(5):490-6

<sup>10</sup> Gürünluoglu et al. 2019: J Burn Care Res. 2019 Jun 21;40(4):444-450

<sup>11</sup> Gürünluoglu et al. 2019: J Burn Care Res. 2019 Apr 26;40(3):302-311


<sup>12</sup> Demircan et al. 2021: Ulus Travma Acil Cerrahi Derg. 2021 Jan;27(1):122-131






## Facts and Figures


### Pain benefits at STSG donor sites



### Pain benefits in cases of burns



## Permeability to water vapor





## Histology

- No inflammatory reaction
- Good vascularized dermal tissue
- Completely formed stratum basale
- Increased epithelization

» 2a° burn: biopsy 14 days after the application of SUPRATHEL®

Uhlig et al.: Burns, 2007 Mar; 33(2):221-9





## The Temporary Second Skin Application



### Wound preparation

Complete removal of all necrotic and nonviable tissue must be performed prior to placement of SUPRATHEL®. The wound bed should be thoroughly cleaned and inspected. All contamination sources should be removed. Second degree burns require thorough debridement to ensure the absence of all contaminated and devitalized tissue.



### Application of SUPRATHEL®

SUPRATHEL® contours well to all parts of the body. After application to the injured site the membrane becomes translucent, allowing for visualization and inspection of the healing process.



### Dressing SUPRATHEL®

SUPRATHEL® should be covered with 1-2 layers of petroleum based gauze which are left intact until the wound is fully healed. Only outer dressings should be changed after inspection of SUPRATHEL® or if they become soiled.



### Removal of SUPRATHEL®

SUPRATHEL® starts to detach from the skin following epithelization and may be removed without causing pain. SUPRATHEL® that is still adhered should be left intact.

## Application with burns on the hands

59 years old, electric burns, mainly superficial partial thickness/deep partial thickness



» Day 1, debridement



» Day 6



» Day 21



» After 8 months

Uhlig et al.: Handchir Mikrochir Plast Chir. 2007 Oct; 39(5):314-9

## Advantages

Due to its excellent adhesion and flexibility it can easily be applied to difficult body parts, such as hands, fingers and toes allowing for early range of motion without disrupting the intimate contact between **SUPRATHEL**<sup>®</sup> and the wound bed.

## Application with large scale burns

38 years old, 95 % BSA, ABSI 13, mainly superficial partial thickness/deep partial thickness



» Day 2, debridement



» Day 2, application of SUPRATHEL®



» Day 7, SUPRATHEL® *in situ*



» 4 weeks after trauma



» 2.5 years after trauma

Uhlig et al.: Osteo trauma care 2007; 15: 2-7

## Treatment concept


- Application of SUPRATHEL® on large areas of burned skin with various depths of burns as primary measure after surgical debridement
- Quick epithelization after a scarless healing of the burns after 8 to 14 days
- Possible transplantation of deeper areas which are not yet healed
- Scarless healing is also possible with large scale burns

## Application on partial thickness burns

36 years old, 90 % BSA, mainly deep partial thickness



» Day 0



» Day 0, debridement and application of SUPRATHEL®



» Day 18

» After 24 months

Kamoliz et al.: Eur Surg (2008) 40/1:19–26

## Advantages

- Quick epithelization of deep partial thickness burns
- SUPRATHEL® serves as primary measure for the application on large scale deep dermal burns
- Available STSG donor sites can be used for obvious full thickness degree burns
- Second transplantation after the specific identification of full thickness burns
- Almost scarless healing even with deep dermal burns

## Application on wounds caused by TEN (Toxic Epidermal Necrolysis)/Lyell's Syndrome

48 years old, TEN (Toxic Epidermal Necrolysis), superficial partial thickness, 80 % BSA



» Day 0



» Day 6



» After 4 weeks

Uhlig, internal study, Marienhospital Stuttgart

### Advantages

- SUPRATHEL® can be applied easily and safely, even to large areas
- Immediate pain relief after application
- Excellent coverage of wounds, no change of SUPRATHEL® needed
- Significantly less effort for nursing staff
- Cost reduction due to high efficiency

## SUPRATHEL<sup>®</sup> packing unit

7.1 x 9.1 in / 18 x 23 cm

7.1 x 3.9 in / 18 x 10 cm

3.5 x 3.9 in / 9 x 10 cm

2.0 x 2.0 in / 5 x 5 cm

Size	Membranes	Order-No.
2.0 x 2.0 in / 5 x 5 cm	5	150505-MDR
3.5 x 3.9 in / 9 x 10 cm	1	110910-MDR
3.5 x 3.9 in / 9 x 10 cm	5	150910-MDR
7.1 x 3.9 in / 18 x 10 cm	1	111810-MDR
7.1 x 3.9 in / 18 x 10 cm	5	151810-MDR
7.1 x 9.1 in / 18 x 23 cm	1	111823-MDR
7.1 x 9.1 in / 18 x 23 cm	5	151823-MDR

Pictures show membranes in their original sizes.



# Scientific publications on SUPRATHEL®

Baartmans MG, Dokter J, den Hollander JC, Kroon AA, Oranje AP. Use of Skin Substitute Dressings in the Treatment of Staphylococcal Scalded Skin Syndrome in Neonates and Young Infants. *Neonatology*. 2011;100(1):9-13.

Blome-Eberwein SA, Amani H, Lozano DD, Gogal C, Boorse D, Pagella P. Burns. A bio-degradable synthetic membrane to treat superficial and deep second degree burn wounds in adults and children – 4 year experience. *Burns*. 2020 Aug 29;50(305-4179(20)30507-6.

Demircan M, Gürnülüoğlu K, Gözükara Bağ HG, Koçbıyık A, Gül M, Üremiş N, Gül S, Gürnülüoğlu S, Türköz Y, Taşçı A. Impaction of the poly-lactic membrane or hydrofiber with silver dressings on the interleukin-6, tumor necrosis factor- $\alpha$ , transforming growth factor- $\beta$ 3 levels in the blood and tissues of pediatric patients with burns. *Ulus Travma Acil Cerrahi Derg*. 2021 Jan;27(1):122-131.

Everett M, Massand S, Davis W, Burkey B, Glat PM. Use of a copolymer dressing on superficial and partial-thickness burns in a paediatric population. *J Wound Care*. 2015 Jul;24(7):S4-8.

Fischer S, Kremer T, Horter J, Schaefer A, Ziegler B, Kneser U, Hirche C. Suprathel® for severe burns in the elderly: Case report and review of the literature. *Burns*. 2016 Aug;42(5):e86-92.

Galati V, Vonthein R, Stang F, Mailaender P, Kisch T. Split thickness skin graft versus application of the temporary skin substitute suprathel in the treatment of deep dermal hand burns: a retrospective cohort study of scar elasticity and perfusion. *Int J Burns Trauma*. 2021 Aug 15;11(4):312-320.

Gürnülüoğlu K, Demircan M, Koç A, Koçbıyık A, Taşçı A, Durmuş K, Gürnülüoğlu S, Gözükara Bağ H. The Effects of Different Burn Dressings on Length of Telomere and Expression of Telomerase in Children With Thermal Burns. *J Burn Care Res*. 2019 Apr 26;40(3):302-311.

Gürnülüoğlu K, Demircan M, Taşçı A, Üremiş MM, Türköz Y, Bağ HG, Akinci A, Bayrakçı E. The effects of two different burn dressings on serum oxidative stress indicators in children with partial burn. *J Burn Care Res*, Jun 21;40(4):444-450.

Harenberg PS, Hrabowski M, Ryssel H, Gazyakan E, Germann G, Engel W, Reichenberger MA. Febrile Ulceronecrotic Mucha-Habermann Disease. *Eplasty*. 2010 Jul 16;10: e53.

Highton L, Wallace C, Shah M. Use of Suprathel® for partial thickness burns in children. *Burns*. 2013 Feb;39(1):136-41.

Hundeshagen G, Collins VN, Wurzer P, Sherman W, Voigt CD, Cambiaso-Daniel J, Nunez Lopez O, Sheaffer J, Herndon DN, Finnerty CC, Branski LK. A Prospective, Randomized, Controlled Trial Comparing the Out-patient Treatment of Pediatric and Adult Partial-Thickness Burns with Suprathel or Mepilex Ag. *J Burn Care Res*. 2018 Feb 20;39(2):261-267.

Kaartinen IS, Kuokkanen HO. Suprathel® causes less bleeding and scarring than Mepilex® Transfer in the treatment of donor sites of split-thickness skin grafts. *J Plast Surg Hand Surg*. 2011 Sep;45(4-5):200-3.

Keck M, Selig HF, Lumenta DB, Kamolz LP, Mittelböck M, Frey M. The use of Suprathel® in deep dermal burns: First results of a prospective study. *Burns*. 2012 May;38(3):388-95.

Lindford AJ, Kaartinen IS, Virolainen S, Vuola J. Comparison of Suprathel® and allograft skin in the treatment of a severe case of toxic epidermal necrolysis. *Burns*. 2011 Nov;37(7):e67-72.

Liodaki E, Schopp BE, Lindert J, Krämer R, Kisch T, Mailänder P, Stang F. Kombination von universellem Antidot und temporärem Hautersatz bei Verätzungen [Combination of a universal antidote and temporary skin substitute for chemical burns: Extended case report]. *Unfallchirurg*. 2015 Sep;118(9):804-7.

Madry R, Struzyna J, Stachura-Kulach A, Drozd L, Bugaj M. Effectiveness of Suprathel® application in partial thickness burns, frostbites and Lyell syndrome treatment. *Pol Przegl Chir*. 2011 Oct 1;83(10):541-8.

Markl P, Prantl L, Schremel S, Babilas P, Landthaler M, Schwarze H. Management of split-thickness donor sites with synthetic wound dressings: results of a comparative clinical study. *Ann Plast Surg*. 2010 Nov;65(5):490-6.

März V, Vogt M. Skin Healing of Deep Second Degree Burn Injuries in Four Individuals Sustained in a Boat Explosion Results after Different Approaches. *Eur. Burn J*. 2020, 1, 191–195.

Merz KM, Sievers R, Reichert B. Suprathel® bei zweitgradig oberflächlichen Verbrennungen im Gesicht [Suprathel® for coverage of superficial dermal burns of the face]. *GMS Verbrennungsmedizin* 2011, Vol 4, ISSN 1869-1412.

Mueller E, Haim M, Petnehazy T, Acham-Roschitz B, Trop M. An innovative local treatment for staphylococcal scalded skin syndrome. *Eur J Clin Microbiol Infect Dis*. 2010 Jul;29(7):893-7.

Nolte SV, Xu W, Rodemann H-P, Rennekampff H-O. Suitability of Biomaterials for Cell Delivery in Vitro. *Osteo trauma care* 2007; 15(1):42-47.

Pfurtscheller K, Trop M. Phototoxic Plant burns: Report of a Case and Review of Topical Wound Treatment in Children. *Pediatr Dermatol*. 2014 Nov-Dec;31(6):e156-9.

Pfurtscheller K, Zobel G, Roedel S, Trop M. Use of Suprathel dressing in a young infant with TEN. *Pediatr Dermatol*. 2008 Sep-Oct;25(5):541-3.

Rapp, M, Uhlig C, Dittel K-K. The Treatment of Mass Burn Casualties Resulting from Mass Disaster. *Osteo trauma care* 2007; 15:8-16.

Rashaan ZM, Krijnen P, Allema JH, Vloemans AF, Schipper IB, Breederfeld RS. Usability and effectiveness of Suprathel® in partial thickness burns in children. *Eur J Trauma Emerg Surg*. 2017 Aug;43(4):549-556.

Ring A, Tilkorn D, Ottomann C, Geomelas M, Steintraesser L, Langer S, Goertz O. Intravital monitoring of microcirculatory and angiogenic response to lactocapromer terpolymer matrix in a wound model. *Int Wound J*. 2011 Apr;8(2):112-7.

Rothenberger J, Constantinescu MA, Held M, Aebersold DM, Stolz A, Tschumi C, Olariu R. Use of a Polylactide-based Copolymer as a Temporary Skin Substitute for a Patient With Moist Desquamation Due to Radiation. *Wounds*. 2016 Jul;28(7):E26-30.

Sari E, Eryılmaz T, Tetik G, Ozakpinar HR, Eker E. Suprathel®-assisted surgical treatment of the hand in a dystrophic epidermolysis bullosa patient. *Int Wound J*. 2014 Oct;11(5):472-5.

Schiefer JL, Rahmanian-Schwarz A, Schaller HE, Manoli T. A Novel Hand-shaped Suprathel simplifies the Treatment of Partial-Thickness Burns. *Adv Skin Wound Care*. 2014 Nov;27(11):513-6.

Schriek K, Ott H, Sinnig M. Paradigm Shift in Treatment Strategies for Second-Degree Burns Using a Caprolactone Dressing (Suprathel®)? A 15-Year Pediatric Burn Center Experience in 2084 Patients. *Eur. Burn J*. 2022, 3, 1–9.

Schwarze H, Küntscher M, Uhlig C, Hierlemann H, Prantl L, Noack N, Hartmann B. Suprathel, a new skin substitute, in the management of donor sites of split-thickness skin grafts: results of a clinical study. *Burns*. 2007 Nov;33(7):850-4.

Schwarze H, Küntscher M, Uhlig C, Hierlemann H, Prantl L, Ottomann C, Hartmann B. Suprathel, a new skin substitute, in the management of partial-thickness burn wounds: results of a clinical study. *Ann Plast Surg*. 2008 Feb;60(2):181-5.

Selig H, Keck M, Lumenta DB, Mittelböck M, Kamolz LP. The use of a polylactide-based copolymer as a temporary skin substitute in deep dermal burns: 1-year follow-up results of a prospective clinical noninferiority trial. *Wound Repair Regen*. 2013 May-Jun;21(3):402-9.

Uhlig C, Rapp M, Dittel KK. Neue Strategien zur Behandlung thermisch geschädigter Hände unter Berücksichtigung des Epithelersatzes Suprathel [New strategies for the treatment of thermally injured hands with regard to the epithelial substitute Suprathel]. *Handchir Mikrochir Plast Chir*. 2007 Oct;39(5):314-9.

Uhlig C, Rapp M, Hartmann B, Hierlemann H, Planck H, Dittel KK. Suprathel-an innovative, resorbable skin substitute for the treatment of burn victims. *Burns*. 2007 Mar;33(2):221-9.

Uhlig C, Hierlemann H, Dittel K-K. Actual Strategies in the Treatment of Severe Burns - Considering Modern Skin Substitutes. *Osteo trauma care* 2007; 15:2-7.



# THE TEMPORARY SECOND SKIN

**PolyMedics Innovations GmbH**

Heerweg 15 D  
73770 Denkendorf | Germany

Phone +49 (0)711 719 500-0  
Fax +49 (0)711 719 500-10  
E-mail [info@polymedics.com](mailto:info@polymedics.com)

[www.polymedics.com](http://www.polymedics.com) | [www.suprathelu.com](http://www.suprathelu.com)

