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Harmonisierte Normen	Kurze Beschreibung
<b>Gemeinsame Spezifikationen</b>	Für das Produkt nicht verfügbar
<b>DIN EN ISO 13485</b>	Medical devices - Quality management systems - Requirements for regulatory purposes
<b>DIN EN 62366-1</b>	Medical devices - Part 1: Application of usability engineering to medical devices
<b>DIN EN ISO 14971</b>	Medical devices – Application of risk management to medical devices
<b>DIN EN ISO 14155</b>	Clinical Investigation of Medical Devices for Human Subjects - Good Clinical Practice
<b>DIN EN ISO 10993-1</b>	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management system
<b>DIN EN ISO 10993-3</b>	Biological evaluation of medical devices - Part 3: Tests for genotoxicity, carcinogenicity and reproductive toxicity
<b>DIN EN ISO 10993-5</b>	Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity
<b>DIN EN ISO 10993-6</b>	Biological evaluation of medical devices – Part 6: Tests for local effects after implantation
<b>DIN EN ISO 10993-10</b>	Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization
<b>DIN EN ISO 10993-11</b>	Biological evaluation of medical devices - Part 11: Tests for systemic toxicity

<b>DIN EN ISO 10993-12</b>	Biological evaluation of medical devices - Part 12: Sample preparation and reference materials
<b>DIN EN ISO 11737-1</b>	Sterilization of medical devices - Requirements for the estimation of population of microorganisms on a product
<b>DIN EN ISO 11737-2</b>	Sterilization of medical devices - Microbiological methods - Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process
<b>DIN EN ISO 11137-1</b>	Sterilization of health care products - Radiation - Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices
<b>DIN EN ISO 11137-2</b>	Sterilization of Health Care Products - Radiation - Part 2: Establishing The Sterilization Dose
<b>DIN EN 556-1</b>	Sterilization of medical devices - Requirements for medical devices to be designated "STERILE" - Part 1: Requirements for terminally sterilized medical devices
<b>DIN EN ISO 11607-1</b>	Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems
<b>DIN EN ISO 11607-2</b>	Packaging for terminally sterilized medical devices - Part 2: Validation requirements for forming, sealing and assembly processes
<b>DIN EN ISO 20417</b>	Medical devices – Information to be supplied by the manufacturer
<b>DIN EN ISO 15223-1</b>	Medical Devices - Symbols To Be Used With Medical Device Labels, Labelling And Information To Be Supplied - Part 1: General Requirements
<b>DIN EN 868-2</b>	Packaging for terminally sterilized medical devices - Part 2: Sterilization wrap - Requirements and test methods
<b>DIN EN 868-5</b>	Packaging for terminally sterilized medical devices - Part 5: Sealable pouches and reels of porous materials and plastic film construction - Requirements and test methods
<b>DIN EN ISO 14698-1</b>	Cleanrooms and associated controlled environments -- Biocontamination control -- Part 1: General principles and methods
<b>DIN EN ISO 14698-2</b>	Cleanrooms and associated controlled environments -- Biocontamination control -- Part 2: Evaluation and interpretation of biocontamination data
<b>ISTA 2a</b>	Partial Simulation Performance Tests - Packaged Products 150 lb (68 kg) or Less
<b>USP &lt;151&gt;</b>	Pyrogen Study
<b>ASTM F1886/F1886M</b>	Standard Test Method for Determining Integrity of Seals for Flexible Packaging by Visual Inspection
<b>ASTM F88/F88M</b>	Standard Test Method for Seal Strength of Flexible Barrier Materials
<b>ASTM F3039</b>	Standard Test Method for Detecting Leaks in Nonporous Packaging or Flexible Barrier Materials by Dye Penetration

## 9. Literaturreferenzen

Baartmans, M. G., Dokter, J., Den Hollander, J. C., Kroon, A. A. & Oranje, A. P. (2011). Use of skin substitute dressings in the treatment of staphylococcal scalded skin syndrome in neonates and young infants. *Neonatology*, 100, 9-13.

Blome-Eberwein, S.A., Amani, H., Lozano, D.D., Gogal, C., Boorse, D., Pagella, P. Burns. A bio-degradable synthetic membrane to treat superficial and deep second degree burn wounds in adults and children – 4 year experience. 2020 Aug 29;S0305-4179(20)30507-6

Demircan, M., Gürnlüoğlu, K., Gözükar, Bag, H.G., Koçbıyık, A., Gül, M., Üremiş, N., Gül, S., Gürnlüoğlu, S., Türköz, Y., Taşçı, A. Impaction of the polylactic membrane or hydrofiber with

silver dressings on the interleukin-6, tumor necrosis factor- $\alpha$ , transforming growth factor-b3 levels in the blood and tissues of pediatric patients with burns. *Ulus Travma Acil Cerrahi Derg.* 2021 Jan;27(1):122-131

Everett, M., Massand, S., Davis, W., Burkey, B. & Glat, P. M. (2015). Use of a copolymer dressing on superficial and partial-thickness burns in a paediatric population. *J Wound Care*, 24, S4-8.

Fischer, S., Kremer, T., Horter, J., Schaefer, A., Ziegler, B., Kneser, U., Hirche, C. Suprathel® for severe burns in the elderly: Case report and review of the literature. *Burns*. 2016 Aug;42(5):e86-92

Gürünlüoğlu, K., Demircan, M., Tasci, A., Uremis, M. M., Turkoz, Y., Bag, H. G., Akinci, A. & Bayrakci, E. (2019). The Effects of Two Different Burn Dressings on Serum Oxidative Stress Indicators in Children with Partial Burn. *J Burn Care Res*, 40, 444-450.

Gürünlüoğlu, K., Demircan, M., Koç, A., Koçbıyık, A., Taşçı, A., Durmuş, K., Gürünlüoğlu, S., Gözükarabağ, H. The Effects of Different Burn Dressings on Length of Telomere and Expression of Telomerase in Children With Thermal Burns. *J Burn Care Res*. 2019 Apr 26;40(3):302-311.

Hakkarainen, T., Koivuniemi, R., Kosonen, M., Escobedo-Lucea, C., Sanz-Garcia, A., Vuola, J., Valtonen, J., Tammela, P., Mäkitie, A., Luukko, K., Yliperttula, M., Kavola, H. *J Control Release*. Nanofibrillar cellulose wound dressing in skin graft donor site treatment. 2016 Dec 28;244(Pt B):292-301.

Harenberg, P. S., Hrabowski, M., Ryssel, H., Gazyakan, E., Germann, G., Engel, H. & Reichenberger, M. A. (2010). CASE REPORT Febrile Ulceronecrotic Mucha-Habermann Disease. *Eplasty*, 10.

Highton, L., Wallace, C., Shah, M. Use of Suprathel® for partial thickness burns in children. *Burns*. 2013 Feb;39(1):136-41

Hundeshagen, G., Collins, V. N., Wurzer, P., Sherman, W., Voigt, C. D., Cambiaso-Daniel, J., Nunez Lopez, O., Sheaffer, J., Herndon, D. N., Finnerty, C. C. & Branski, L. K. (2018). A Prospective, Randomized, Controlled Trial Comparing the Outpatient Treatment of Pediatric and Adult Partial-Thickness Burns with Suprathel or Mepilex Ag. *J Burn Care Res*, 39, 261-267

Kaartinen, I. S. & Kuokkanen, H. O. (2011). Suprathel® causes less bleeding and scarring than Mepilex® Transfer in the treatment of donor sites of split-thickness skin grafts. *J Plast Surg Hand Surg*, 45, 200-3.

Keck, M., Selig, H. F., Lumenta, D. B., Kamolz, L. P., Mittlbock, M. & Frey, M. (2012). The use of Suprathel® in deep dermal burns: first results of a prospective study. *Burns*, 38, 388-95.

Lindford, A. J., Kaartinen, I. S., Virolainen, S. & Vuola, J. (2011). Comparison of Suprathel® and allograft skin in the treatment of a severe case of toxic epidermal necrolysis. *Burns*, 37, e67-72.

- Liodaki, E., Schopp, B.E., Lindert, J., Krämer, R., Kisch, T., Mailänder, P., Stang, F. Kombination von universellem Antidot und temporärem Hautersatz bei Verätzungen [Combination of a universal antidote and temporary skin substitute for chemical burns: Extended case report]. Unfallchirurg. 2015 Sep;118(9):804-7.
- Madry, R., Struzyna, J., Stachura-Kulach, A., Drozd, L. & Bugaj, M. (2011). Effectiveness of Suprathel(R) application in partial thickness burns, frostbites and Lyell syndrome treatment. Pol Przegl Chir, 83, 541-8.
- Markl, P., Prantl, L., Schreml, S., Babilas, P., Landthaler, M. & Schwarze, H. (2010). Management of split-thickness donor sites with synthetic wound dressings: results of a comparative clinical study. Ann Plast Surg, 65, 490-6.
- Merz, K. M., Sievers, R., Reichert, B. (2011). Suprathel® for coverage of superficial dermal burns of the face. GMS Verbrennungsmedizin, 4.
- Mueller, E., Haim, M., Petnehazy, T., Acham-Roschitz, B. & Trop, M. (2010). An innovative local treatment for staphylococcal scalded skin syndrome. Eur J Clin Microbiol Infect Dis, 29, 893-7.
- Pfurtscheller, K., Zobel, G., Roedl, S. & Trop, M. (2008). Use of Suprathel dressing in a young infant with TEN. Pediatr Dermatol, 25, 541-3.
- Pfurtscheller, K. & Trop, M. (2014). Phototoxic plant burns: report of a case and review of topical wound treatment in children. Pediatr Dermatol, 31, e156-9.
- Rahmanian-Schwarz, A., Beiderwieden, A., Willkomm, L.M., Amr, A., Schaller, H.E., Lotter, O. A clinical evaluation of Biobrane(®) and Suprathel(®) in acute burns and reconstructive surgery. Burns. 2011 Dec;37(8):1343-8
- Rashaan, Z. M., Krijnen, P., Allema, J. H., Vloemans, A. F., Schipper, I. B. & Breederveld, R. S. (2017). Usability and effectiveness of Suprathel((R)) in partial thickness burns in children. Eur J Trauma Emerg Surg, 43, 549-556.
- Rothenberger, J., Constantinescu, M. A., Held, M., Aebersold, D. M., Stolz, A., Tschumi, C. & Olariu, R. (2016). Use of a Polylactide-based Copolymer as a Temporary Skin Substitute for a Patient With Moist Desquamation Due to Radiation. Wounds, 28, E26-30.
- Sari, E., Eryilmaz, T., Tetik, G., Ozakpinar, H. R. & Eker, E. (2014). Suprathel((R)) -assisted surgical treatment of the hand in a dystrophic epidermolysis bullosa patient. Int Wound J, 11, 472-5.
- Schiefer, J.L., Rahmanian-Schwarz, A., Schaller, H.E., Manoli, T. A Novel Hand-shaped Suprathel simplifies the Treatment of Partial-Thickness Burns. Adv Skin Wound Care. 2014 Nov;27(11):513-6

- Schwarze, H., Kuntscher, M., Uhlig, C., Hierlemann, H., Prantl, L., Noack, N. & Hartmann, B. (2007). Suprathel, a new skin substitute, in the management of donor sites of split-thickness skin grafts: results of a clinical study. *Burns*, 33, 850-4.
- Schwarze, H., Kuntscher, M., Uhlig, C., Hierlemann, H., Prantl, L., Ottomann, C. & Hartmann, B. (2008). Suprathel, a new skin substitute, in the management of partial-thickness burn wounds: results of a clinical study. *Ann Plast Surg*, 60, 181-5.
- Selig, H. F., Keck, M., Lumenta, D. B., Mittlbock, M. & Kamolz, L. P. (2013). The use of a polylactide-based copolymer as a temporary skin substitute in deep dermal burns: 1-year follow-up results of a prospective clinical noninferiority trial. *Wound Repair Regen*, 21, 402-9.
- Uhlig, C., Hierlemann, H., Dittel, K.-K. (2007). Actual Strategies in the Treatment of Severe Burns - Considering Modern Skin Substitutes. *Osteo trauma care* 15, 2-7.
- Uhlig, C., Rapp, M. & Dittel, K. K. (2007a). [New strategies for the treatment of thermally injured hands with regard to the epithelial substitute Suprathel]. *Handchir Mikrochir Plast Chir*, 39, 314-9.
- Uhlig, C., Rapp, M., Hartmann, B., Hierlemann, H., Planck, H. & Dittel, K. K. (2007b). Suprathel-an innovative, resorbable skin substitute for the treatment of burn victims. *Burns*, 33, 221-9.

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