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Normas armonizadas	Breve descripción
Common specifications	Non available for the product
DIN EN ISO 13485	Medical devices - Quality management systems - Requirements for regulatory purposes
DIN EN 62366-1	Medical devices - Part 1: Application of usability engineering to medical devices
DIN EN ISO 14971	Medical devices – Application of risk management to medical devices
DIN EN ISO 14155	Clinical Investigation of Medical Devices for Human Subjects - Good Clinical Practice
DIN EN ISO 10993-1	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management system
DIN EN ISO 10993-3	Biological evaluation of medical devices - Part 3: Tests for genotoxicity, carcinogenicity and reproductive toxicity
DIN EN ISO 10993-5	Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity
DIN EN ISO 10993-6	Biological evaluation of medical devices – Part 6: Tests for local effects after implantation
DIN EN ISO 10993-10	Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization
DIN EN ISO 10993-11	Biological evaluation of medical devices - Part 11: Tests for systemic toxicity
DIN EN ISO 10993-12	Biological evaluation of medical devices - Part 12: Sample preparation and reference materials
DIN EN ISO 11737-1	Sterilization of medical devices - Requirements for the estimation of population of microorganisms on a product
DIN EN ISO 11737-2	Sterilization of medical devices - Microbiological methods - Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process
DIN EN ISO 11137-1	Sterilization of health care products - Radiation - Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices
DIN EN ISO 11137-2	Sterilization of Health Care Products - Radiation - Part 2: Establishing The Sterilization Dose
DIN EN 556-1	Sterilization of medical devices - Requirements for medical devices to be designated "STERILE" - Part 1: Requirements for terminally sterilized medical devices
DIN EN ISO 11607-1	Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems
DIN EN ISO 11607-2	Packaging for terminally sterilized medical devices - Part 2: Validation requirements for forming, sealing and assembly processes
DIN EN ISO 20417	Medical devices – Information to be supplied by the manufacturer
DIN EN ISO 15223-1	Medical Devices - Symbols To Be Used With Medical Device Labels, Labelling And Information To Be Supplied - Part 1: General Requirements
DIN EN 868-2	Packaging for terminally sterilized medical devices - Part 2: Sterilization wrap - Requirements and test methods
DIN EN 868-5	Packaging for terminally sterilized medical devices - Part 5: Sealable pouches and reels of porous materials and plastic film construction - Requirements and test methods

DIN EN ISO 14698-1	Cleanrooms and associated controlled environments -- Biocontamination control -- Part 1: General principles and methods
DIN EN ISO 14698-2	Cleanrooms and associated controlled environments -- Biocontamination control -- Part 2: Evaluation and interpretation of biocontamination data
ISTA 2a	Partial Simulation Performance Tests - Packaged Products 150 lb (68 kg) or Less
USP <151>	Pyrogen Study
ASTM F1886/F1886M	Standard Test Method for Determining Integrity of Seals for Flexible Packaging by Visual Inspection
ASTM F88/F88M	Standard Test Method for Seal Strength of Flexible Barrier Materials
ASTM F3039	Standard Test Method for Detecting Leaks in Nonporous Packaging or Flexible Barrier Materials by Dye Penetration

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Baartmans, M. G., Dokter, J., Den Hollander, J. C., Kroon, A. A. & Oranje, A. P. (2011). Use of skin substitute dressings in the treatment of staphylococcal scalded skin syndrome in neonates and young infants. *Neonatology*, 100, 9-13.

Blome-Eberwein, S.A., Amani, H., Lozano, D.D., Gogal, C., Boorse, D., Pagella, P. Burns. A bio-degradable synthetic membrane to treat superficial and deep second degree burn wounds in adults and children – 4 year experience. 2020 Aug 29;S0305-4179(20)30507-6

Demircan, M., Gürünlüoğlu, K., Gözükar, Bag, H.G., Koçbıyık, A., Gül, M., Üremiş, N., Gül, S., Gürünlüoğlu, S., Türköz, Y., Taşçı, A. Impaction of the polylactic membrane or hydrofiber with silver dressings on the interleukin-6, tumor necrosis factor- α , transforming growth factor-b3 levels in the blood and tissues of pediatric patients with burns. *Ulus Travma Acil Cerrahi Derg.* 2021 Jan;27(1):122-131

Everett, M., Massand, S., Davis, W., Burkey, B. & Glat, P. M. (2015). Use of a copolymer dressing on superficial and partial-thickness burns in a paediatric population. *J Wound Care*, 24, S4-8.

Fischer, S., Kremer, T., Horter, J., Schaefer, A., Ziegler, B., Kneser, U., Hirche, C. Suprathel® for severe burns in the elderly: Case report and review of the literature. *Burns*. 2016 Aug;42(5):e86-92

Gürünlüoğlu, K., Demircan, M., Tasci, A., Uremis, M. M., Turkoz, Y., Bag, H. G., Akinci, A. & Bayrakci, E. (2019). The Effects of Two Different Burn Dressings on Serum Oxidative Stress Indicators in Children with Partial Burn. *J Burn Care Res*, 40, 444-450.

Gürünlüoğlu, K., Demircan, M., Koç, A., Koçbıyık, A., Taşçı, A., Durmuş, K., Gürünlüoğlu, S., Gözükar Bağ, H. The Effects of Different Burn Dressings on Length of Telomere and Expression of Telomerase in Children With Thermal Burns. *J Burn Care Res*. 2019 Apr 26;40(3):302-311.

- Hakkarainen, T., Koivuniemi, R., Kosonen, M., Escobedo-Lucea, C., Sanz-Garcia, A., Vuola, J., Valtonen, J., Tammela, P., Mäkitie, A., Luukko, K., Yliperttula, M., Kavola, H. J Control Release. Nanofibrillar cellulose wound dressing in skin graft donor site treatment. 2016 Dec 28;244(Pt B):292-301.
- Harenberg, P. S., Hrabowski, M., Ryssel, H., Gazyakan, E., Germann, G., Engel, H. & Reichenberger, M. A. (2010). CASE REPORT Febrile Ulceronecrotic Mucha-Habermann Disease. Eplasty, 10.
- Highton, L., Wallace, C., Shah, M. Use of Suprathel® for partial thickness burns in children. Burns. 2013 Feb;39(1):136-41
- Hundeshagen, G., Collins, V. N., Wurzer, P., Sherman, W., Voigt, C. D., Cambiaso-Daniel, J., Nunez Lopez, O., Sheaffer, J., Herndon, D. N., Finnerty, C. C. & Branski, L. K. (2018). A Prospective, Randomized, Controlled Trial Comparing the Outpatient Treatment of Pediatric and Adult Partial-Thickness Burns with Suprathel or Mepilex Ag. J Burn Care Res, 39, 261-267
- Kaartinen, I. S. & Kuokkanen, H. O. (2011). Suprathel((R)) causes less bleeding and scarring than Mepilex((R)) Transfer in the treatment of donor sites of split-thickness skin grafts. J Plast Surg Hand Surg, 45, 200-3.
- Keck, M., Selig, H. F., Lumenta, D. B., Kamolz, L. P., Mittlbock, M. & Frey, M. (2012). The use of Suprathel((R)) in deep dermal burns: first results of a prospective study. Burns, 38, 388-95.
- Lindford, A. J., Kaartinen, I. S., Virolainen, S. & Vuola, J. (2011). Comparison of Suprathel(R) and allograft skin in the treatment of a severe case of toxic epidermal necrolysis. Burns, 37, e67-72.
- Liodaki, E., Schopp, B.E., Lindert, J., Krämer, R., Kisch, T., Mailänder, P., Stang, F. Kombination von universellem Antidot und temporärem Hautersatz bei Verätzungen [Combination of a universal antidote and temporary skin substitute for chemical burns: Extended case report]. Unfallchirurg. 2015 Sep;118(9):804-7.
- Madry, R., Struzyna, J., Stachura-Kulach, A., Drozd, L. & Bugaj, M. (2011). Effectiveness of Suprathel(R) application in partial thickness burns, frostbites and Lyell syndrome treatment. Pol Przegl Chir, 83, 541-8.
- Markl, P., Prantl, L., Schreml, S., Babilas, P., Landthaler, M. & Schwarze, H. (2010). Management of split-thickness donor sites with synthetic wound dressings: results of a comparative clinical study. Ann Plast Surg, 65, 490-6.
- Merz, K. M., Sievers, R., Reichert, B. (2011). Suprathel® for coverage of superficial dermal burns of the face. GMS Verbrennungsmedizin, 4.
- Mueller, E., Haim, M., Petnehazy, T., Acham-Roschitz, B. & Trop, M. (2010). An innovative local treatment for staphylococcal scalded skin syndrome. Eur J Clin Microbiol Infect Dis, 29, 893-7.

- Pfurtscheller, K., Zobel, G., Roedl, S. & Trop, M. (2008). Use of Suprathel dressing in a young infant with TEN. *Pediatr Dermatol*, 25, 541-3.
- Pfurtscheller, K. & Trop, M. (2014). Phototoxic plant burns: report of a case and review of topical wound treatment in children. *Pediatr Dermatol*, 31, e156-9.
- Rahmanian-Schwarz, A., Beiderwieden, A., Willkomm, L.M., Amr, A., Schaller, H.E., Lotter, O. A clinical evaluation of Biobrane® and Suprathel® in acute burns and reconstructive surgery. *Burns*. 2011 Dec;37(8):1343-8
- Rashaan, Z. M., Krijnen, P., Allema, J. H., Vloemans, A. F., Schipper, I. B. & Breederveld, R. S. (2017). Usability and effectiveness of Suprathel® in partial thickness burns in children. *Eur J Trauma Emerg Surg*, 43, 549-556.
- Rothenberger, J., Constantinescu, M. A., Held, M., Aebersold, D. M., Stolz, A., Tschumi, C. & Olariu, R. (2016). Use of a Polylactide-based Copolymer as a Temporary Skin Substitute for a Patient With Moist Desquamation Due to Radiation. *Wounds*, 28, E26-30.
- Sari, E., Eryilmaz, T., Tetik, G., Ozakpinar, H. R. & Eker, E. (2014). Suprathel® -assisted surgical treatment of the hand in a dystrophic epidermolysis bullosa patient. *Int Wound J*, 11, 472-5.
- Schiefer, J.L., Rahmanian-Schwarz, A., Schaller, H.E., Manoli, T. A Novel Hand-shaped Suprathel simplifies the Treatment of Partial-Thickness Burns. *Adv Skin Wound Care*. 2014 Nov;27(11):513-6
- Schwarze, H., Kuntscher, M., Uhlig, C., Hierlemann, H., Prantl, L., Noack, N. & Hartmann, B. (2007). Suprathel, a new skin substitute, in the management of donor sites of split-thickness skin grafts: results of a clinical study. *Burns*, 33, 850-4.
- Schwarze, H., Kuntscher, M., Uhlig, C., Hierlemann, H., Prantl, L., Ottomann, C. & Hartmann, B. (2008). Suprathel, a new skin substitute, in the management of partial-thickness burn wounds: results of a clinical study. *Ann Plast Surg*, 60, 181-5.
- Selig, H. F., Keck, M., Lumenta, D. B., Mittlbock, M. & Kamolz, L. P. (2013). The use of a polylactide-based copolymer as a temporary skin substitute in deep dermal burns: 1-year follow-up results of a prospective clinical noninferiority trial. *Wound Repair Regen*, 21, 402-9.
- Uhlig, C., Hierlemann, H., Dittel, K.-K. (2007). Actual Strategies in the Treatment of Severe Burns - Considering Modern Skin Substitutes. *Osteo trauma care* 15, 2-7.
- Uhlig, C., Rapp, M. & Dittel, K. K. (2007a). [New strategies for the treatment of thermally injured hands with regard to the epithelial substitute Suprathel]. *Handchir Mikrochir Plast Chir*, 39, 314-9.
- Uhlig, C., Rapp, M., Hartmann, B., Hierlemann, H., Planck, H. & Dittel, K. K. (2007b). Suprathel-an innovative, resorbable skin substitute for the treatment of burn victims. *Burns*, 33, 221-9.

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