Autologous Cell Harvesting Device

# Instructions for use

ENGLISH

FRANÇAIS

ESPAÑOL

PORTUGUÉS

DEUTSCH

ITALIANO

NEDERLANDS

TÜRK

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**avita** medical

CE  
1639

Autologous Cell Harvesting Device

# Instructions for use

ReCell® should be used only by licensed health care professionals trained in the use of ReCell®

**Warning:**

The ReCell Autologous Cell Harvesting Device is internally powered by four non-replaceable AA batteries (1.5v). The device should not be used in the presence of flammable anaesthetic mixtures. Do not incinerate on disposal. The performance of the device may be affected by sources of electromagnetic radiation and if any malfunctions are noted, all possible sources of electromagnetic radiation must be removed before further use.

# TABLE OF CONTENTS

<b>A</b>	<b>BACKGROUND</b>	<b>3</b>
A1	DESCRIPTION .....	3
A2	INTENDED USE / INDICATIONS.....	3
A3	CONTRAINDICATIONS .....	3
A4	WARNINGS .....	3
A5	PRECAUTIONS .....	3
A6	ADVERSE REACTIONS.....	3
A7	MEANING OF SYMBOLS.....	4
A8	DOSAGE .....	4
A9	HOW SUPPLIED .....	4
A10	STORAGE .....	4
A11	DISPOSAL.....	4
<b>B</b>	<b>RECELL TREATMENT</b>	<b>5</b>
B1	MATERIALS.....	5
B2	RECELL DEVICE SET UP.....	5
B3	WOUND BED PREPARATION .....	6
B4	STEP-BY-STEP INSTRUCTIONS FOR PREPARING THE REGENERATIVE EPITHELIAL SUSPENSION (RES).....	6
<b>C</b>	<b>AFTERCARE</b>	<b>10</b>
C1	INITIAL DRESSING .....	10
C2	SUBSEQUENT DRESSINGS .....	10
C3	AFTERCARE PRECAUTIONS.....	10
C4	SCAR MANAGEMENT.....	10
<b>D</b>	<b>TROUBLESHOOTING</b>	<b>11</b>

# A BACKGROUND

## A1 DESCRIPTION

ReCell® is a single-use, stand-alone, battery-operated, autologous cell harvesting device. It contains a proprietary enzymatic cell disaggregation and activation solution, as well as delivery solutions, sterile surgical instruments, and actuators. The ReCell device enables thin split-thickness skin samples to be processed to produce a Regenerative Epithelial Suspension™ (RES™) for immediate delivery onto a prepared wound surface.

The RES contains a mixed population of cells obtained from the disaggregation of the skin sample including primarily keratinocytes and fibroblasts, but also melanocytes, Langerhans cells, and epidermal basal cells. The preservation of melanocytes is important for restoring natural pigmentation to the recipient area.

The Enzyme used to process the cells is a biological agent and as such may have slight variations in colour and texture.

## A2 INTENDED USE / INDICATIONS

ReCell is intended to be used to disaggregate cells from a patient's split-thickness skin sample and to collect these cells for reintroduction to the patient. The cells can be used for autologous application to the prepared wound bed as determined by the physician such as for the treatment of burns, or other acute wounds.

## A3 CONTRAINDICATIONS

- ReCell is contraindicated for patients with wounds that are clinically infected or necrotic.
- ReCell should not be used to prepare Regenerative Epithelial Suspension for application to patients with a known hypersensitivity to trypsin or compound sodium lactate solution.
- The skin sample collection procedure specified for use of ReCell should not be used with patients having a known hypersensitivity to anaesthetics, adrenaline/epinephrine, povidine-iodine, or chlorhexidine solutions.

## A4 WARNINGS

- The Regenerative Epithelial Suspension produced with ReCell should only be applied to the patient from whom the original skin sample was taken (autologous use only).
- ReCell is provided to the healthcare professional sterile and is intended for single use. Do not reuse, freeze or re-sterilise device components.
- Do not use ReCell or device components if packaging is damaged or there are signs of tampering.
- Do not use ReCell or device components if the date of use is beyond the stated expiration date on the packaging.
- ReCell components should be handled using aseptic technique.
- If a skin sample is harvested and processed according to these instructions, it should only require between 15 and 30 minutes of contact with the Enzyme. Contact in excess of 60 minutes is not recommended.
- Contaminated materials and waste must be disposed of using appropriate biohazard waste receptacles.
- The separation Enzyme is derived from animal tissue and, although strict controls have been implemented in the manufacturing process to minimize the risk of pathogen contamination, a small risk of contamination exists and absolute freedom from infectious agents cannot be guaranteed.
- The ReCell processing unit is internally powered by four non-replaceable AA batteries. The device should not be used in the presence of flammable materials and must not be incinerated on disposal.

## A5 PRECAUTIONS

- Protective eyewear and other protective clothing should be worn.
- For optimum cell viability, the skin sample should be processed immediately after harvesting.
- The ReCell device is for single use only. Do not reuse, freeze or re-sterilise any items within the device.
- Do not use the device if there is evidence of container tampering or damage.

## A6 ADVERSE REACTIONS

Any adverse reaction or suspected adverse reaction related to ReCell should immediately be reported to Avita Medical.

## A7 MEANING OF SYMBOLS

The packaging system is labelled with various symbols. These symbols are internationally harmonised and define certain characteristics of the product and the manufacturing process:

This symbol states that the product is for single use only

This symbol states that the user should refer to the accompanying instructions for use

This symbol states that the date adjacent is the date of manufacture

This symbol states that the date adjacent is the expiry date of the product

This symbol states the manufacturer of the product

This symbol states that the temperature adjacent specifies the upper limit of storage temperature

This symbol states that the product or components within have been sterilised using ethylene oxide

This symbol states that the product or components within have been sterilised using gamma irradiation

This symbol states that the product or components within have been sterilised using steam

## A8 DOSAGE

ReCell is supplied as a single use device. The contents of each device are sufficient to prepare up to 24 ml of Regenerative Epithelial Suspension which can be used to cover an acute wound area of up to and including 1920 cm<sup>2</sup>.

## A9 HOW SUPPLIED

The ReCell Device consists of:

- 1 x Processing Unit with built-in heating mechanism
  - 1 x removable sterile tray
  - 1 x removable cell strainer
- 1 x sealed vial of Enzyme
- 1 x 10-ml vial of sterile water
- 4 x 10-ml vial of Buffer
- 7 x 10-ml syringes
- 2 x sharp needles
- 4 x blunt drawing-up needles
- 1 x additional cell strainer
- 2 x disposable surgical scalpel
- 1 x packet sterile syringe labels
- 4 x spray nozzles

Component sterilisation and testing

- The Processing Unit and needles have been sterilised by ethylene oxide.
- The Enzyme has undergone filtration and terminal sterilisation by gamma irradiation.
- The scalpel and spray nozzles have been sterilised by gamma irradiation.
- The syringes have been sterilised by either ethylene oxide or gamma irradiation.
- The Buffer and sterile water have been sterilised using steam.

## A10 STORAGE

Upon receiving ReCell, examine the packaging for external signs of damage. If the external packaging or the packaging for any of the individual components appears damaged, contact your local representative immediately. Do not use any components of the device if the packaging appears damaged. If returning ReCell, ensure all original packaging and components are returned with the device.

ReCell, including the Enzyme, may be stored at controlled room temperature, not to exceed 27 °C. Testing indicates that the Enzyme remains active for at least 21 months at this temperature. The Enzyme is labelled with an expiry date of 21 months after manufacture.

Do not open or use ReCell outside of the expiration date listed on the packaging.

## A11 DISPOSAL

- ReCell and all individual components are intended for single use. ReCell components are not reusable and should be discarded after single use. Reuse may lead to infection or disease transmission.
- Follow local regulations for proper disposal.
- Contaminated materials and waste must be disposed of using appropriate biohazard receptacles.
- ReCell contains batteries and electrical components – DO NOT INCINERATE
- A procedure for removal of processing unit battery/electronics is as follows:
  - Take proper Biohazard precautions when handling the used processing unit.
  - Remove the processing unit top cover. Set top cover aside.
  - Remove inner tray and set aside.
  - Open inner main tray by pressing both sides of the outer housing simultaneously.
  - Verify that the parts are separated (inner main tray and outer housing). If the parts of the inner tray and outer housing are not separated, a small, flat-blade screwdriver may be used to assist in releasing the inner and outer parts.
  - Lift the inner tray to expose battery compartment.
  - Remove the batteries and the electronics and dispose of them in the appropriate waste streams.
  - Dispose of the remaining components in accordance with the appropriate methods.

## B RECELL TREATMENT

### B1 MATERIALS

The following materials and instruments will be needed during the ReCell procedure.

- Surgical gloves and a suitable sterile drape
- Protective eyewear and garments
- Skin preparation (antiseptic) solution
- Local anaesthetic with adrenaline where not contraindicated
- Appropriate wound dressings. See "Aftercare" below for details
- 1 or 2 x fine-point (long nosed) forceps of choice
- Skin harvesting instrument of choice, e.g., Zimmer® Dermatome (Zimmer Orthopedic Surgical Products, Inc., USA), Silver's knife, Humby knife, Dermablade® (Personna® American Safety Razor Co, USA)
- Wound bed preparation tool of choice
- Clock or timer to monitor incubation time

### B2 RECELL DEVICE SET UP

Select and prepare sterile and non-sterile work areas. Note that the ReCell device contains both sterile and non-sterile components.

Using standard aseptic technique, set up a sterile surgical field.

- Remove the Processing Unit from the sterile packaging and place it in the sterile field.
- Open the Processing Unit and note the removable inner white plastic insert. This insert acts as a sterile tray for use in preparing and scraping the skin sample.

### PERFORM SELF-TEST

Perform the self-test to verify the processing unit is functioning correctly.

- Test the processing unit to ensure functionality by pressing the button marked (?). All lights should illuminate during the self-test. When the unit has completed the self-test (this takes approximately 30 seconds), it will beep once and the green 'ready' light (✓) will illuminate to indicate that the Processing Unit is functioning correctly. If lights do not illuminate, or the red light (!) illuminates, do not use the device. The unit will automatically turn off after 1 minute if Enzyme heating is not initiated.
- If the device turns off after self-test, additional self-tests may be run.
- Do not press the run button (▶) at this time.

### A - PREPARE ENZYME (COMPONENT SET A)


- In the non-sterile work area, remove the cover from the vial marked Enzyme to expose the injection diaphragm. Wipe the diaphragm with a sterile alcohol wipe and allow to dry (optional).
- Connect a sharp, sterile needle to a sterile 10-ml syringe and draw up the entire volume of sterile water.
- Inject the entire volume of sterile water into the Enzyme vial. DO NOT USE Buffer at this stage as this may inhibit the Enzyme action.
- Tip the vial gently to mix the solution until dissolved. Do not shake and use care to avoid foaming.
- Draw the entire volume of dissolved Enzyme back into the syringe.
- Using aseptic technique, dispense the entire volume of Enzyme into the left-hand well of the Processing Unit (Well A). Discard syringe and needle.

### B - PREPARE BUFFER SOLUTION (COMPONENT SET B)

Introduce the following items into the sterile field.

- 2 x 10-ml syringes
- 1 x 19ga sharp needle
- 1 x sheet sterile syringe labels
- 1 x cell strainer
- 2 x disposable surgical scalpel

The four (4) Buffer vials are to remain outside the sterile field

- Apply the "BUFFER" label to one of the new 10-ml syringes. The "BUFFER" syringe will be used several times to draw Buffer from the vials. Set aside within the sterile field.
- Apply the "UNFILTERED SUSPENSION" label to the other 10-ml syringe. The "UNFILTERED SUSPENSION" syringe will be used several times to collect unfiltered cell suspension from the tray and dispense into the cell strainer. Set aside within the sterile field.
- It is important that these syringes are used only for their intended, labelled purpose and that they remain sterile.
- Remove the cover from a vial marked Buffer. Wipe the diaphragm of the vial with a sterile alcohol wipe and allow it to dry (optional).
- Attach the sharp, sterile needle to "BUFFER" syringe and draw up the entire volume of Buffer (approximately 10 ml) from the vial.
- Dispense the entire volume of Buffer into the empty centre well of the Processing Unit (Well B). The entire 10 ml of Buffer will be used to rinse skin samples.
- Set the "BUFFER" syringe and needle aside within the sterile field for later use.

**C - PREPARE DELIVERY ITEMS (COMPONENT SET C)**

Introduce all items into the sterile field.

- 4 x 10-ml syringes
- 4 x blunt drawing-up needles
- 4 x spray nozzles

**B3 WOUND BED PREPARATION**

- Clean, vascularized wound bed – To optimise the treatment, the Regenerative Epithelial Suspension should only be applied to a clean, vascularised wound bed with no remaining necrotic tissue. This can be achieved with either dermabrasion using a rotating diamond-head burr, laser ablation, sharp dissection, or other alternative techniques depending on the nature of the wound.
- Infection free – The Regenerative Epithelial Suspension must not be used in the presence of any contamination or infection, as initial re-epithelialisation and long term viability are highly dependent on the absence of infection. Prophylactic antibiotics may be prescribed if the patient is at risk of contamination or infection. Wound swabs for up-to-date microbiology are recommended 48 hours prior to the planned surgery.
- Pinpoint bleeding – The wound bed should be prepared so that dermis is exposed and pinpoint bleeding is observed. This may be accomplished via several methods, e.g., ablative laser, fractional laser, or mechanical abrasion. As much viable dermis as possible should be preserved. If tissue damage is present (e.g., burns) accurate debridement to the level of viable tissue is essential; all necrotic tissue must be removed.
- When ReCell is used for treatment of acute, full-thickness wounds, combination with split-thickness skin graft and/or dermal replacement technology is recommended.

**B4 STEP-BY-STEP INSTRUCTIONS FOR PREPARING THE REGENERATIVE EPITHELIAL SUSPENSION (RES)**

**1. Take Skin sample**

**Skin sample Type**

It is essential that the skin sample harvested is a thin, split-thickness skin sample that penetrates to the dermis and leaves pinpoint bleeding at the donor site. The thickness of the skin sample will vary with the body site and patient age and should be in the range of 0.15 to 0.20 mm (0.006 to 0.008 in). The use of a dermatome, Silver’s knife, Humby knife or DermaBlade® is recommended.

**Size of Skin Sample**

Choose the appropriate skin sample size for the application. Each square centimetre of skin sample area can create up to 1 ml of RES for treatment of an area of up to 80 cm<sup>2</sup>. Each 6 cm<sup>2</sup> (3 cm x 2 cm) skin sample can yield approximately 6 ml of RES; each ReCell device can process up to 4, 6 cm<sup>2</sup> skin samples for a maximum of 24 ml of RES. This can be used to treat an area of approximately 1,920 cm<sup>2</sup>. Sample guidance for skin sample size is shown in the following table.

Treatment Area (Including donor site)	Skin Sample Size
up to 80 cm <sup>2</sup>	1 cm x 1 cm (1 cm <sup>2</sup> )
up to 320 cm <sup>2</sup>	2 cm x 2 cm (4 cm <sup>2</sup> )
up to 480 cm <sup>2</sup>	3 cm x 2 cm (6 cm <sup>2</sup> )
up to 960 cm <sup>2</sup>	2 ea. 3 cm x 2 cm (12 cm <sup>2</sup> )
up to 1440 cm <sup>2</sup>	3 ea. 3 cm x 2 cm (18 cm <sup>2</sup> )
up to 1920 cm <sup>2</sup>	4 ea. 3 cm x 2 cm (24 cm <sup>2</sup> )

**Choice of Donor Site**

It is essential the donor site is clean, of appropriate depth, and shows no evidence of surrounding cellulitis or infection. Choose a donor site of glabrous tissue when creating suspension for glabrous tissue regeneration.

**Harvesting the skin sample**

Using the preferred instrument such as a dermatome, Silver’s knife, or Humby knife, take a split-thickness shave biopsy of the donor site of thickness 0.15 to 0.20 mm (0.006 to 0.008 in). Use the table above to estimate the skin sample size needed, or calculate by taking 1/80 of the total treatment area (wound plus donor sites areas).

The following settings are recommended for obtaining a skin sample of the proper thickness using a Zimmer dermatome.

	Children	Adult
Blade guard (width)	2.5 cm	4.5 cm
Dermatome Setting	6	8

Clean the donor site with antiseptic solution such as povidine-iodine or chlorhexidine. Allow the antiseptic to dry before removing with sterile saline (antiseptic solutions may be cytotoxic and as such, may affect cell viability if left on the skin sample site).

If desired, infiltrate the subcutaneous tissue with a tumescent solution of choice, to provide a firmer surface and anaesthesia for taking the skin sample. Ensure that anaesthetic is not injected intradermally.

The donor site area may be lubricated, for instance with sterile mineral oil, to ease travel of the dermatome.

The donor site may be treated similarly to the primary wound. Once haemostasis has been achieved, a small amount of RES may be applied to the donor site wound, and dressed as described in the Aftercare section below.

Due to the thick keratin layer found on glabrous skin, it is necessary to take two shaves over the same site in these areas. Discard the first sample and process the second skin sample to create the cell suspension.

## 2. Heat Enzyme

Verify that the Enzyme has been transferred to Well A. The processing unit will quickly overheat if the run button (▶) is pressed before the Enzyme has been placed in the well. Any malfunctioning of the unit, including overheating, will be indicated by the red light (!) illuminating. Should this occur, use another ReCell device and contact your local representative to arrange the return or replacement of the unit.

Press the run button (▶) to heat the Enzyme in Well A. If the device is ready (✓) then heating will commence. If more than one minute has passed since the last self-test, a self-test will automatically run, followed immediately by heating of the well. The orange warming light will illuminate when warming begins and the Enzyme will be heated and maintained at the approximately 37 °C.

## 3. Incubate the Skin Samples

When the orange warming light turns off and the green (✓) illuminates the Enzyme has reached its target temperature. This will take approximately 3 minutes. The orange light will flash from time to time, indicating that the heating element has been activated to maintain temperature.

Place 1 or 2 skin samples into the heated Enzyme for 15 to 20 minutes to allow extra-cellular matrix breakdown. Skin samples (6 cm<sup>2</sup>) may be incubated 2 at a time. Up to 4 samples may be processed using a single device, 2 initially, then 2 following those. Keep the other skin samples moist in sterile gauze moistened with sterile saline. If the skin samples are thick, they may require longer incubation. Each sample may be incubated for up to 60 minutes however the Processing Unit will maintain the enzyme temperature for only 75 minutes.

## 4. Draw up Buffer

This step may be performed whilst the skin sample is incubating.

Remove the cover from a vial marked Buffer. Wipe the diaphragm of the vial with a sterile alcohol wipe and allow it to dry (optional). Using aseptic technique and the "BUFFER" syringe with sharp needle, draw up the required volume of Buffer from the Buffer vial. Use one millilitre of Buffer per square centimetre of the first skin sample and add 0.5 ml Buffer to allow for loss during processing. The following table provides example surface areas to be treated, skin sample sizes needed, volumes of Buffer to use, and approximate resultant RES volumes.

Place the syringe with Buffer in the sterile field for use in steps 8 and 9, below.

Surface Area to be Treated per Syringe	Skin Sample Size Needed	Starting Volume of Buffer	Approximate Resultant RES Volume
Up to 80 cm <sup>2</sup>	1 cm <sup>2</sup> (1 cm x 1 cm)	1.5 ml	1.0 ml
Up to 160 cm <sup>2</sup>	2 cm <sup>2</sup> (1 cm x 2 cm)	2.5 ml	2.0 ml
Up to 320 cm <sup>2</sup>	4 cm <sup>2</sup> (2 cm x 2 cm)	4.5 ml	4.0 ml
Up to 480 cm <sup>2</sup>	6 cm <sup>2</sup> (3 cm x 2 cm)	6.5 ml	6.0 ml

ReCell RES may be used alone or in combination with meshed, split-thickness grafts, or dermal substitutes. Each 1 ml of RES may be used to treat of up to 80 square centimetres when used on acute partial-thickness wounds.

## 5. Test for Cell Disaggregation

After 15 to 20 minutes, remove the skin sample(s) from the heated enzyme solution with sterile forceps and place the skin sample dermal side down on the sterile tray. Gently scrape the epidermis with the scalpel to test if cells disaggregate, i.e., epidermal cells easily come off. If the cells do not come off freely, return the skin sample(s) to the heated Enzyme for a further 5 to 10 minutes and then repeat the test scrape. When the cells scrape off freely, proceed to the next step.

After approximately 60 minutes, an alarm will sound and will sound each minute for 15 minutes. At 75 minutes, the Processing Unit will turn off and stop heating the enzyme. Incubation of a skin sample for more than 60 minutes is not recommended.

## 6. Rinse Skin Samples

Upon a successful test scrape, briefly rinse the skin sample in the middle well (Well B) containing the Buffer to rinse off and deactivate the residual Enzyme. Return the skin sample to the sterile tray.

## 7. Incubate Additional Skin Samples

If processing 3 or 4 samples, place these into Well A for 15-20 minutes.

### 8. Scrape Cells from the Skin sample

With the skin sample dermal side down on the sterile tray, apply a few drops of Buffer from the previously filled "BUFFER" syringe onto the skin sample. Using the forceps to anchor the skin sample, gently scrape the epidermal surface with the blade of the scalpel and collect in the pool of buffer. Once the epidermis has been scraped away, scrape the remaining dermis more rigorously. Continue scraping until the dermis has nearly disintegrated.

### 9. Rinse and Aspirate; Draw up unfiltered suspension

Use the remaining Buffer in the "BUFFER" syringe to rinse the scalpel and tray, collecting the unfiltered suspension into one corner of the tray. Tilt the tray as necessary. Set the "BUFFER" syringe aside for later use. Using the "UNFILTERED SUSPENSION" syringe, collect and draw up the unfiltered suspension. Using the drawn-up suspension, rinse the tray. Draw up and rinse several times to maximize cell collection. Finally, draw up the suspension into the syringe.

### 10. Filter Cells

Dispense the unfiltered cell suspension into the cell strainer in Well C. Set the "UNFILTERED SUSPENSION" syringe aside, within the sterile field, for use with unfiltered suspensions from the remaining skin samples.

### 11. Draw up Regenerative Epithelial Suspension (RES)

Attach a blunt needle to a new 10-ml syringe. Carefully remove the cell strainer, tapping the cell strainer over the well to release any residual drops of RES. Draw up the filtered RES from Well C. There is a conical point in the centre of the bottom of Well C to aid in drawing up all of the RES. Return the cell strainer to Well C.

## 12. Process Sample 2

Process the second skin sample in the same manner as the first. Using the "BUFFER" syringe and sharp needle, draw up the appropriate volume of Buffer for the skin sample size from the Buffer vials (Step 4). Repeat steps 8-11 (scrape, rinse, collect, filter, and draw up) for the second skin sample to produce a suspension. Use the "UNFILTERED SUSPENSION" syringe to collect, rinse, and draw up the unfiltered suspension from the tray and transfer to the cell strainer. Use a new 10-ml syringe and blunt needle to draw up the RES. After sample 2 has been processed, replace the cell strainer in Well C with the second cell strainer provided.

## 13. Process Samples 3 and 4

After the samples 3 and 4 have incubated in the enzyme for 15-20 minutes remove them from Well A. Perform the test scrape (step 5) then rinse if ready (step 6). Process first one, then the other by repeating steps 8-11 for each sample. Use the "BUFFER" syringe to draw up fresh Buffer and the "UNFILTERED SUSPENSION" to collect, rinse, and draw up the suspension from the tray and transfer to the cell strainer. Use new 10-ml syringes to draw up the RES from Well C for each sample. Return the cell strainer to Well C.

### 14. Apply RES to Wound Bed

Prior to applying the RES, ensure the dressings are cut and prepared for immediate application. The primary dressing should be fixed or held at the lower aspect of the wound prior to applying the RES. Section C, Aftercare, provides information on dressing selection and use.

The RES can be sprayed using the spray nozzles provided, dripped onto the wound or introduced under the primary dressing using a blunt drawing-up needle.

The choice of application method depends on the volume of RES and the size and location of the wound(s) (refer to table below). Some RES may be reserved for application to the donor site.

RES Volume per Syringe	Surface Area to be Treated	Example Wound Site	Recommended Application Method
1.0 ml	Up to 80 cm <sup>2</sup>	Palms	Drip or beneath dressing
2.0 ml	Up to 160 cm <sup>2</sup>	Face, Neck	Spray or Drip
6.0 ml	Up to 480 cm <sup>2</sup>	Trunk, Extremities	Spray

The minimum volume of RES required for spray application is approximately 2 ml.

#### Spray Application

Remove the needle from the syringe containing the RES. Attach the spray nozzle supplied to the syringe using firm pressure. Invert the syringe several times prior to the application to ensure an even suspension. Check that the aperture of the attached spray nozzle faces the wound. Hold the spray applicator approximately 10 cm from the most elevated point of the wound and in a position such that the first drop of RES falls onto the wound surface. Apply moderate pressure to the plunger of the syringe. Start spraying at the most elevated part of the wound so that any run-off helps to cover the more dependent areas of the wound. A fine mist of RES should be delivered to the wound surface. To cover a larger area, carefully move the spray applicator in one continuous motion from one side of the wound to the other as you spray.

#### Drip Application

Do not remove the blunt needle from the syringe containing the RES. Invert the syringe several times prior to application to ensure an even suspension. Starting at the most elevated point of the wound, carefully drip the RES onto the wound surface.

#### Application under Primary Dressing

If introducing the RES under a dressing do not remove the blunt needle from the syringe containing the RES. Invert the syringe several times prior to application to ensure an even suspension. Place the cut dressing over the wound and gently introduce the needle under the dressing and introduce the RES. Larger wounds may require introducing the needle and RES at several points to ensure complete coverage.

Note: The fibrin in the prepared wound bed provides an ideal environment to cell adhesion. Many, but not all, of the delivered cells will adhere to the wound. It is normal for some of the RES to run off the wound. A well-prepared suspension has enough cells to treat the wound, allowing for run off.

## C AFTERCARE

The following information, precautions, and notes provide guidelines for care after ReCell. Discuss appropriate aftercare with your Avita representative and provide patient with the "After ReCell" guidance brochure.

### C1 INITIAL DRESSING

After applying the RES, cover the wound with a non-adherent, non-absorbent, small pore dressing such as Telfa™ Clear (Covidien, USA) wound dressing or similar. In some cases it may be more practical to apply the primary dressing first and introduce the RES under the dressing using a blunt needle. Always follow the instructions as set by the dressing manufacturer. Dry dressings such as Telfa™ Clear wound dressing may be applied moist at the direction of the healthcare professional by lightly soaking the dressing in sterile saline before dressing the wound. The dressing may be fixed to the wound with surgical glue, sutures, or staples, as necessary.

Place a secondary dressing over the primary dressing to provide protection to the wound and absorb any exudate. This should be a moderately absorbent, minimally adherent, low shear, readily removable dressing. One example is a paraffin gauze layer placed over the primary dressing, covered with a moist saline or Povidine-Iodine compress and a crepe bandage.

### C2 SUBSEQUENT DRESSINGS

The outer dressings and compression may need to be changed if exudate levels are high, however the primary dressing should remain in place for 6-8 days, or as clinically indicated. Take care to protect the primary dressing during secondary dressing changes. The primary dressing will loosen and lift as new epidermis is formed and should not be removed from areas to which it is still adhered. IT IS ESSENTIAL THAT PRIMARY DRESSING REMOVAL IS ATRAUMATIC. ANY DRESSING NOT EASILY REMOVED SHOULD BE SOAKED WITH AN AQUEOUS OR OIL-BASED SOLUTION PRIOR TO REMOVAL TO PREVENT TRAUMA. Once the primary dressing has been removed, an appropriate protective dressing such as Jelonet® (Smith & Nephew, UK) or Mepitel® (Mölnlycke, Sweden) should be applied to protect the wound surface.

Do not use dry dressings as protection over an area of punctate blistering, as dried exudate could cause newly regenerated skin to adhere to the dressing, leading to potential injury upon dressing removal. Instead use (for instance) a greasy or paraffin gauze dressing until any blistering or open areas resolve.

Any signs or symptoms of infection or impaired healing at this stage should be recorded and addressed.

### C3 AFTERCARE PRECAUTIONS

- Patients should take necessary precautions to prevent the treated area from getting wet while the wound is still open.
- Do not disrupt the primary dressing for a minimum of 6-8 days. Ensure that primary dressing removal is atraumatic - Do not forcibly remove the primary dressing.
- Patients should be advised to avoid trauma to the wound and dressing, including off-loading as appropriate. Up to two additional weeks may be needed after initial closure of the treated area for the newly regenerated skin to mature and become robust. During this time protective dressings must be worn, particularly on extremities.
- Use of known cytotoxic medication (for instance, silver sulfadiazine) is contraindicated for areas treated using ReCell.
- Patients and caregivers should be provided with adequate information and materials for appropriate protection against re-injury during healing and maturation of the treated area.
- Patients should be advised to refrain from strenuous activity.
- Patients should avoid direct sun exposure for at least four weeks following treatment

### C4 SCAR MANAGEMENT

When the wound has healed, the patient should be advised to continue to protect the area from any surface trauma and to avoid direct sun for at least four weeks. Regular use of sun block and twice-daily massage with a non-oily skin moisturizer is recommended.

The patient should be advised that the wound area will change over the subsequent weeks and months. The pigmentation and skin texture will continue to mature and improve during this time and the final result may take up to 12 months to be achieved.

Follow-up procedures should follow standard protocols for the specific injury and treatment given.

## D TROUBLESHOOTING

### **Enzyme powder does not dissolve completely**

Make sure that the Enzyme is mixed well with the sterile water by inverting the vial several times. Often a small amount of particulate matter remains undissolved in the reconstituted solution. This does not reduce the activity of the Enzyme.

Do not use Buffer to dissolve the Enzyme as it may interfere with the Enzyme action.

### **Skin sample is too large, too thick or too thin**

Take particular care when harvesting the skin sample. It should be a thin (0.15 to 0.20 mm) split-thickness shave biopsy with just a very thin section of dermis (see previous instructions for dermatome settings). The skin sample of the appropriate thickness will ensure successful disaggregation of cells. The maximum size of skin sample recommended for use with the ReCell device is 3 cm by 2 cm.

If the skin sample is too large (greater than the maximum recommended), cut it into a smaller size and discard the excess.

If the skin sample is too thick, cut the skin sample into 1 cm by 1 cm pieces before placing in the heated Enzyme. If the cells cannot be disaggregated, repeatedly return the skin sample to the heated Enzyme for a further 5 to 10 minutes, up to a maximum of 60 minutes of total time. If the cells still do not scrape off freely it may be necessary to take another thin, split-thickness skin sample from a DIFFERENT donor site and repeat the process using a new ReCell device.

If the skin sample is too thin, you should take another skin sample from a DIFFERENT donor site and repeat the process.

### **Buffer added to Enzyme vial**

If Buffer is mistakenly added to the Enzyme vial, instead of sterile water, the Enzyme activity may be inhibited. If Buffer is mixed with the Enzyme powder, the Enzyme should be discarded and a new ReCell device used.

### **Difficult Cell Disaggregation**

Ensure that the heating element is switched on. The green light (✓) will illuminate when the ReCell device is switched on and ready for use. The orange warming light will illuminate when the device is warming. Disaggregation of the cells will take longer if the skin sample is too large or thick. See above for advice.

### **Nozzle blocked**

If the RES is not easily sprayed or does not come out at all, the nozzle attached to the syringe may be blocked. Use one of the other nozzles provided or consider dripping the suspension on to the wound.

### **Insufficient treatment area coverage**

If RES is lost in the application process and sufficient coverage of the treatment area was not achieved, take another skin sample and repeat the process with a new ReCell device to create additional RES and complete the treatment.

For further information regarding ReCell Autologous Cell Harvesting Device, contact your local sales representative at:

**AVITA Medical Americas, LLC**  
**(Legal Manufacturer)**  
**28159 Avenue Stanford, Suite #200**  
**Valencia, CA 91355-2271, USA**  
**Tel: +1 661 367 9170**  
**Fax: +1 661 367 9180**  
**Email: support.am@avitamedical.com**

**Emergo Europe**  
**Westervoortsedijk 60**  
**6827 AT Arnhem**  
**The Netherlands**  
**Tel: +31 0 70 345 8570**  
**Fax: +31 0 70 346 7299**

The ReCell® Autologous Cell Harvesting Device is subject to issued and pending patent and design applications of Avita Medical.  
US Patent Nos. 9,029,140 & 9,078,741; European Patent Nos. 1 357 922 & 2 343 079 AT, BE, DE, ES, FR, GB, IT, NL, PT, SE, TR;  
Japan Patent No. 5214085; Hong Kong Patent No. 1057713

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© Avita Medical 2015